



OneLABz

User Manual (0.1)

By

TechXect Limited

www.techxect.com

Table of drafts

Login Screen

Viewing Test Results Without Login

Application Admin (Login with administrator user)

Admin

Add Lab

Labs

Add home sample

Home samples

Custom Notification

Admin settings

Currency

Measuring unit

Items

Test settings

Tests

Test Packages

Categories

Parameter type

Collection type

Specimen

Parameter result template

Sample status

Discount option

Parameter Interpretation

Admin parameter

Lab settings

Collection center

Employee

Doctor

Access control

Add roles

Roles

Settings

Add collection centre

Collection centres

Add chemical

Chemicals listing

Lab Configurations

Patient

Add Patient

Patients

Stock

Manage Stock

Stock Listing

Quality Control System

Add Quality Control values

Quality controls values

Sample Collection

Collect Sample

Adding Test

Payment details

Samples

Test Result

Add test result

Test result

All test results

Pending perform

Pending approval

Expense

Add Expense

Expenses

Expense Report

Add expense type

Expense types

Reports

Patient dues amount report

Doctor shares Report

Daily summary report

Doctor Share Summary Report

Test summary report

Payment detail report

Company rates report

Payment received report


Stock consumption report

Discount option report

OneLABz onboarding


Login Screen

Enter your given username and password into the appropriate boxes, then click the 'Login' button to gain access to the OneLABz application.



Viewing Test Results Without Login

If a patient does not have an account or does not wish to log in, test results can be viewed by clicking the **View test result** link located on the Sign-in screen. They can then enter the required information and click the 'View' button to see their test results.



An example of a report print and the result is as below

[Print](#)



EcoMed
ADVANCED MEDICAL DIAGNOSTICS

EcoMed Lab®

EcoMed (SL) Ltd., 19 Off Upper Kingharman Road, (100 meters from Ministry of Tourism), Brookfields, Freetown, Sierra Leone.
 Email: info@ecomeds.l Website: http://onelabz.techco.net/
 Phone/WhatsApp: +232 90 999 529



Lab id: 50934 Collection centre: Main lab Sample collected on: 22-Jan-2020 11:29 Sample collected by: Fazeela Kanwal Performed on: 22-Jan-2020 11:36 Performed by: Hafiz Abdul Hafeez	Patient name: Zafar Iqbal Father/husband name: Gulzar ali Age / Sex: 28(Y)/Male Marital status: 03448794849 Phone number: sammad.abbas@btsol.com Email: 145/S-L Address: Respected Dr CNIC: Consultant:	
--	---	--

Serology				Printed On: 30-Mar-2020 16:00:42
Tests	Parameter	Result	Reference value	Unit
Blood Grouping	(Blood Grouping)*	A-VE(Negative)	---	

Dr. Farah Naz Tahir Amin
 Consultant Pathologist:

Electronically verified report, no signature(s) required.
 Lab values should be interpreted by a physician in the context of clinical picture.


Jr. Farah naz Tahir amin MBBS, M. PHIL., Consultant Pathologist	Hafiz Abdul Hafeez Medical Lab Tech. N.I.H Islamabad	Shahryar Habib ullah MBBS
--	---	------------------------------

Note: All the tests are performed on the most advanced, highly sophisticated, appropriate, and state of the art instruments with highly sensitive chemicals under strict conditions and with all care and diligence. However, the above results are NOT the DIAGNOSIS and should be correlated with clinical findings, patient's history, signs and symptoms and other diagnostic tests. Lab to lab variation may occur. This document is NEVER challengeable at any PLACE/COURT and in any CONDITION.

Powered by TechMest

Application Admin (Login with administrator user)

In order to set up the application for lab use, the system administrator should add certain details in the software e.g. Lab employees, their login details, tests performed, collection centers, chemicals, their stock, etc.




Dashboard

No. of Samples Per Day

In this graph daily collected samples are displayed. Data is displaying the current month.


Dashboard



Filter on graph

Admin can view graphs against all Collection centers, so we provided Collection center filters and date filters


Dashboard



to search desired data.

Top tests performed in 30 days

A bar graph is displaying the performed test. Data is displaying 30 days.




Admin can view graphs against all Collection centers, so we provided Collection center filters and date filters.

Lowest stock items


Stock updates are displayed in this graph.

We are providing filters to stock against the different labs.




No of created patient

Created patients are displayed in this graph.




Using filters can check patients against different collection centers.



Top patients by outstanding dues


Patient monthly outstanding dues are displayed in this graph.




The filter is provided to search patient dues against any collection center.

Top patients by revenue

Monthly patient revenue is displayed in this graph.




Can filter revenue against desired collection center.




Top patients by tests

Top patients are displayed in this graph who did more tests in a month. Date filters and collection center filters are provided to patients.



Top doctors by share in last 30 days


Doctor share is displayed in this graph.




Check the share against any doctor using a filter.

Top doctors by referrals

When collecting samples and referring to any doctor that is displayed in this graph.



Using filters can view how many times doctors are referred.




Collection center today expense and revenue

In this graph, the Collection center can see its daily expenses and revenue.




This graph can filter against its collection center.




Collection center expense and revenue per day

In this graph, users can view records of one month.




Can filter against the collection center and month.




Top referred employees revenue

Admin can view the referral employee's revenue in this graph.



Admin can view employee's revenue using filters against employee and date.



MR Monthly target

In this graph, admin can view the target achieved and the monthly target of the marketing representative.



Can use filters to view employees and dates.

MR Monthly target

Filters

Lab: Inova lab

From date: 08/05/2021

To date: 08/06/2021


- Mr Huanan
- Mr-ahsan Ali
- Sania Hussan
- Employee Testing
- MR Bilal
- MR_Hamid Hamid
- MR_Testing Testing
- MR_Faisa Jamal

<--Please select a value-->


Close Apply

MR tests by category

This graph display the category count used in a month.




Users can use filters to filter against categories and dates.




MR tests by collection center

Admin can view the count of tests performed in a month against any collection center.




Top Category Type in Month

In this graph view those categories which are more used in a month.




Can filter which collection center is using categories.




Payments and receipt for the collection center

In this graph, admin can view the total received and remaining amounts.




Can view the received amount for a specific collection center using filters.




Collection centers revenue in the month

Admin can view the monthly revenue of any collection center.




Can filter collection center against date.




Top test category in month

The comparative graph in which admin can view three months data.




Filter used to view monthly against any collection center and of any date.




MR Performance in month


The comparative graph in which can compare three months performance of a marketing representative employee.



Can filter against employee and date.



Admin



Test Settings

Hover the mouse on **Admin** click on **Test settings**, all are present in the tabular form. When clicking

The Test Settings page includes the following elements:

- Navigation:** Dashboard, Sample collection, Test result, Patient, Accounting, Stock, Quality control, Reports, Settings, Access control, Admin.
- Test Settings Section:** Tests, Test packages, Test categories, Test parameter type, Specimen collection type, Test specimen, Test parameter result template.
- Export Options:** Copy, CSV, Excel, PDF, Print.
- Search:** Search by [Please select a value-->] Search value [] Search [] Reset []
- Table:**

Test id	Code	Category	Name	Specimen	Description	Charges (Rs.)	Cre	Me
28696	AARBCC	Special Chemistry	Testing Test	Blood		5000.000	16-0	16-0
977	172	Serology	CRP (Qualitative)	Blood		400.000	18-0	18-0
976	171	Special Chemistry	TPHA Elisa	Blood		1400.000	09-01-2020	Techxect Lab 09-0
975	169	Clinical Chemistry	MAGNESIUM	Blood		300.000	09-01-2020	Techxect Lab 09-0
974	168	Clinical Chemistry	Urine Sugar	Urine		50.000	07-01-2020	Techxect Lab 07-0
973	159	Serology	RA Factor Qualitative	Blood		400.000	01-01-2020	Techxect Lab 01-0
972	158	Clinical Chemistry	Blood glucose Report F...	Blood		100.000	01-01-2020	Techxect Lab 01-0
971	157	Molecular Biology	Viral Marker(Hepatitis ...	Blood		350.000	29-12-2019	Techxect Lab 29-1
970	156	Molecular Biology	Urine For Ovulation Test	Urine		600.000	03-12-2019	Techxect Lab 03-1
969	155	Chemical Pathology	Semen Analysis Report	Semen		300.000	03-12-2019	Techxect Lab 03-1
968	154	Microbiology	Urine Complete Exami...	Urine		150.000	03-12-2019	Techxect Lab 03-1
967	151	Special Chemistry	Prolactin Level	Blood		1000.000	26-11-2019	Techxect Lab 26-1
966	150	Special Chemistry	Progesterone Level	Blood		1000.000	26-11-2019	Techxect Lab 26-1
- Dropdown Menu:** Add doctor, Doctors, Add employee, Employees, Add test, Tests, Add currency, Currencies, Add measuring unit, Measuring units, Admin settings, Test settings, Lab settings.

on **Test settings** by default it is on the **Tests** screen.

Add new Test

On the left side, a button is displaying and clicking on it to add a new Test. Tests that are performed in your laboratory, should be added to the system before they can be conducted for a patient.

There are three sections on the 'Add test' page i.e. Test (details), Test parameters, and Test chemicals/tools (used to perform the test).

To add a test, use the 'Test' menu and select the 'Add test' option. Fill in the form with the required test's details.

Report templates:

We are providing three types of report templates.

- General report
- PCR report
- Custom report

You can use whichever you want.

ONELABS Dashboard Sample collection Test result Patient Accounting Stock Quality control Reports Settings Access control Admin Inova lab Inova lab

Add test

Test

Lab * Inova lab

Test code * Test code

ETA in day(s) - hours(s) - minutes(s) * 0 0 0

Test category * <--Please select a value-->

Test charges (Rs.) * 0.000

Test name * Test name

Test specimen * Test specimen

Test report template * <--Please select a value-->

Manage stock at sample collection

Required approval for print

Notes

Method

Back to tests listing

Manage stock at sample level :

If you want to manage stock at sample collection please check the checkbox **Manage stock at sample collection**.

ONELABS Dashboard Sample collection Test result Patient Accounting Stock Quality control Reports Settings Access control Admin Inova lab Inova lab

Add test

Test

Lab * Inova lab

Test code * Test code

ETA in day(s) - hours(s) - minutes(s) * 0 0 0

Test category * <--Please select a value-->

Test charges (Rs.) * 0.000

Test name * Test name

Test specimen * Test specimen

Test report template * <--Please select a value-->

Manage stock at sample collection

Required approval for print

Notes

Method

Back to tests listing

Test parameters [Add test parameter](#)

Parameter name	Male				Female				Child				Child2				Infant			
	Critical min	Min	Max	Critical max	Critical min	Min	Max	Critical max	Critical min	Min	Max	Critical max	Critical min	Min	Max	Critical max	Critical min	Min	Max	Critical max

Test chemicals / tools [Test chemicals / tools](#)

Activate Windows Go to Settings to activate Windows.

Note: If you want to manage stock when the test is performed then please uncheck it.

The screenshot shows the 'Add test' form in the ONELABz system. The form includes fields for Lab, Test code, ETA, Test report template, Test description, Comments, Test category, Test charges, Test name, Test specimen, and Required approval for print. The 'Manage stock at sample collection' checkbox is highlighted with a red box and is currently checked.

Approval for print:

If you want to print without the approval of the test then uncheck the checkbox **required approval for print**.

The screenshot shows the 'Add test' form in the ONELABz system. The 'Required approval for print' checkbox is highlighted with a red box and is currently unchecked.

If you want to print with approval of the test then check the checkbox **required approval for print**

The screenshot shows the 'Add test' form in the ONELABz system. The 'Required approval for print' checkbox is highlighted with a red box and is currently checked.

Test parameters should be added in the 'Test parameters' section. Use the 'Add test parameter' button to add a new parameter record for the test. To remove a parameter from the test, use the 'Remove' button at the end of the parameter record.

Chemicals or tools consumed for a test are added in the 'Test chemicals/tools' section. Use the 'Add test chemical' button to add a new chemical/tool record for the test. To remove a chemical from the test, use the 'Remove' button at the end of the chemical record. Save the test details using the 'Save' button.

Tests

All tests which are added in OneLABz can be viewed on this screen.


Test id	Code	Category	Name	Specimen	Description	Charges (Rs.)
816	172	Serology	CRP (Qualitative)	Blood		400.000
815	171	Special Chemist...	TPHA Elisa	Blood		1400.000
814	169	Clinical Chemis...	MAGNESIUM	Blood		300.000
813	168	Clinical Chemis...	Urine Sugar	Urine	Urine Sugar	50.000
812	159	Serology	RA Factor Qualitative	Blood		400.000
811	158	Clinical Chemis...	Blood glucose Report Fastl...	Blood	Fasting & Random Sugar	100.000
810	157	Molecular Biolo...	Viral Marker(Hepatitis B&C)	Blood		350.000
809	156	Molecular Biolo...	Urine For Ovulation Test	Urine	Urine For Ovulation Test	600.000
808	155	Chemical Patho...	Semen Analysis Report	Semen	Semen Analysis Report	300.000
807	154	Microbiology	Urine Complete Examination	Urine	Urine Complete Examination	1500.000
806	151	Soecial Chemist...	Prolactin Level	Blood	Prolactin Level	1000.000

Add new Test package:

If the lab wants to give a package against more than one test, then click on the Add **New Test package** button.


There are two sections on the 'Add test package' page i.e. package(details) and add tests that are in the package. To add a test package, use the 'Test setting' menu and select the 'Test settings' and click on the 'Test package' option. Fill all required information and test's package details.

Save the test package details using the 'Save' button.




Test Packages:

All Test packages which are added in OneLABz can be viewed on this screen.



Add test category:

There is one section on the 'Add test category' page i.e. category(details).



To add a test category, use the 'Admin' menu and select the 'Test settings' and click on the 'Test category' option. Fill in the form with the test's category details.

Save the test category details using the 'Save' button.

Test categories:

All Test categories which are added in OneLABz can be viewed on this screen.

The screenshot shows the 'Test settings' page in OneLABz. The 'Test categories' tab is highlighted in red. Below the tabs, there are export options (Copy, CSV, Excel, PDF, Print) and a search bar. A table lists various test categories with columns for Test category id, Test category, Test category description, Created date, Created by, Modified date, Modified by, and Lab.

Test category id	Test category	Test category description	Created date	Created by	Modified date	Modified by	Lab
2661	Covid19	Covid19	18-03-2021	Inova lab Inova lab	18-03-2021	Inova lab Inova lab	Mahmood
2660	Special Chemistry	Special Chemistry	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
2659	None	None	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
2658	MRI	MRI	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
2657	XRAY	XRAY	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
2656	Chemical Pathology	Chemical Pathology	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
2655	Serology	Serology	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
2654	Toxicology	Toxicology	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
2653	Endocrinology	Endocrinology	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
2652	Bacteriology	Bacteriology	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
2651	Parasitology	Parasitology	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
2650	Histopathology	Histopathology	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
2649	Molecular Biology	Molecular Biology	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
2648	Heamatology	Heamatology	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
2647	Clinical Chemistry	Clinical Chemistry	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
2646	Microbiology	Microbiology	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
2645	Sugar	Sugar	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
2644	DNA	DNA	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
2643	Blood Test	Blood Test	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood

Add Test parameter type:

There is one section on the 'Add test parameter type' page i.e. Test parameter (details).To

The screenshot shows the 'Add test parameter type' form. It has a dropdown menu for 'Lab' (set to 'Inova lab'), a text input field for 'Test parameter type name', and another text input field for 'Test parameter type description'. There is a 'Save' button at the bottom left and a 'Back to parameter types listing' link at the top right.

add a test parameter, use the 'Admin' menu and select the 'Test settings' and click on the 'Test parameter type' option. Fill in the form with the test's parameter details.

Save the test parameter details using the 'Save' button.

Test parameter type:


All Test parameter types which are added in OneLABz can be viewed on this screen.

The screenshot shows the 'Test settings' page in OneLABz. The 'Test parameter type' tab is highlighted in red. Below the tabs, there are export options (Copy, CSV, Excel, PDF, Print) and a search bar. A table lists various test parameter types with columns for Parameter type id, Parameter type, Parameter type description, Created date, Created by, Modified date, Modified by, and Lab.

Parameter type id	Parameter type	Parameter type description	Created date	Created by	Modified date	Modified by	Lab
623	None	None	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
622	BioChemical	BioChemical	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
621	Physical	Physical	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
620	Chemical	Chemical	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
599	None	None	20-01-2021	Super Admin	20-01-2021	Super Admin	Child_lab_inova
598	BioChemical	BioChemical	20-01-2021	Super Admin	20-01-2021	Super Admin	Child_lab_inova

Add specimen collection type:


There is one section on the 'Add specimen collection type' page i.e. specimen collection type (details).



To add a specimen collection type, use the 'Admin' menu and select the 'Test settings' and click on the 'Specimen collection type' option. Fill in the form with the test's specimen details. Save the specimen collection type details using the 'Save' button.

Specimen collection type:

All Test specimen collection types which are added in OneLABz can be viewed on this screen.




Specimen collection type id	Specimen collection type	Specimen collection type description	Created date	Created by	Modified date	Modified by	Lab
338	Brought to lab	Brought to lab	21-04-2021	Super Admin	21-04-2021	Super Admin	Mahmood
337	Taken in lab	Taken in lab	21-04-2021	Super Admin	21-04-2021	Super Admin	Mahmood
326	Brought to lab	Brought to lab	20-04-2021	Super Admin	20-04-2021	Super Admin	Child_lab_inova
325	Taken in lab	Taken in lab	20-04-2021	Super Admin	20-04-2021	Super Admin	Child_lab_inova
247	Brought to lab	Brought to lab	03-03-2021	Super Admin	03-03-2021	Super Admin	Inova lab
246	Taken in lab	Taken in lab	03-03-2021	Super Admin	03-03-2021	Super Admin	Inova lab

Add test specimen:


To add a specimen collection type, use the 'Admin' menu and select the 'Test settings' and click on the 'Test specimen' option. Fill in the form with the test's specimen details.

Save the test specimen details using the 'Save' button.



Test Specimen:

All Test specimens which are added in OneLABz can be viewed on this screen.




The screenshot shows the 'Test settings' page in the OneLABz application. The 'Test specimen' tab is selected and highlighted with a red box. Below the navigation tabs, there are export options (Copy, CSV, Excel, PDF, Print) and a search bar. A table lists 12 test specimens with columns for Specimen id, Specimen name, Specimen description, Created date, Created by, Modified date, Modified by, and Lab. The table shows various specimen types like None, Oxygen, Urine, and Blood, all created by 'Super Admin' on '20-01-2021'. At the bottom, it indicates 'Showing 25 records per page. 1-12 Of 12' and a pagination control.

Specimen id	Specimen name	Specimen description	Created date	Created by	Modified date	Modified by	Lab
651	None	None	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
650	Oxygen	Oxygen	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
649	Urine	Urine	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
648	Blood	Blood	20-01-2021	Super Admin	20-01-2021	Super Admin	Mahmood
624	None	None	20-01-2021	Super Admin	20-01-2021	Super Admin	Child_lab_inova
623	Oxygen	Oxygen	20-01-2021	Super Admin	20-01-2021	Super Admin	Child_lab_inova
622	Urine	Urine	20-01-2021	Super Admin	20-01-2021	Super Admin	Child_lab_inova
621	Blood	Blood	20-01-2021	Super Admin	20-01-2021	Super Admin	Child_lab_inova
476	None	None	20-01-2021	Super Admin	20-01-2021	Super Admin	Inova lab
475	Oxygen	Oxygen	20-01-2021	Super Admin	20-01-2021	Super Admin	Inova lab
474	Urine	Urine	20-01-2021	Super Admin	20-01-2021	Super Admin	Inova lab
473	Blood	Blood	20-01-2021	Super Admin	20-01-2021	Super Admin	Inova lab

Add test parameter template:

To add a test parameter template, use the 'Admin' menu and select the 'Test settings' and click on the 'Test parameter template option. Fill in the form with the test's parameter template details.


Save the test parameter template details using the 'Save' button.



The screenshot shows the 'Add test parameter template' form in the OneLABz application. The form has two main sections: 'Test parameter result template' and 'Test parameter result value'. The first section includes a dropdown for 'Lab' (set to 'Inova lab'), a text input for 'Template name', and a checkbox for 'For comment'. The second section has a text input for 'Parameter result' and an 'Add result value' button. A 'Save' button is located at the bottom left of the form. A 'Back to test parameter template listing' link is in the top right corner.

Test parameter templates:


All Test parameters templates which are added in OneLABz can be viewed on this screen.



Sample status:


To add a sample status, use the 'Admin' menu and select the 'Test settings' and click on the 'Sample status' option. Fill in the form with the sample template details. The select the color that the lab wants to use as an indicator.

Save the sample status details using the 'Save' button.



Sample statuses:

All Test parameters templates which are added in OneLABz can be viewed on this screen.



Discount Option

Admin adds a discount option from Test settings for use on different screens.

Discount option

Lab *
Al nasar Lab

Discount option name *
Discount option name

Discount option description
Discription

Discount Percentage
0

Save

[Back to discount option listing](#)

Discount Option listing

Discount option listing displaying in this screen user can copy, print, download in CSV, PDF, and excel form. Searching and sorting are also implemented here.

Test settings

Export
Copy CSV Excel PDF Print

Search by <--Please select a value--> Search value Search value Search Reset

Discount option id	Discount option name	Discount option description	Discount%	Created date	Created by	Modified date	Modifie
53	Discount Option	Discount Option	10	06-07-2021	Al nasar Lab Al nasar L...	06-07-2021	Al nasar L

Showing 25 records per page. 1-1 Of 1

<<Start < 1 > End>>

Lab Settings

There are three types of lab settings, the screenshot is attached below.

Lab settings

Collection center Employee Doctor

Export
Copy CSV Excel PDF Print

Search by <--Please select a v... Search value Search value Search Reset

Collection centre id	Code	Collection center	Address	Description	Max disc per test (%)	Max disc per d	Notification settings	Head office	Enabled	Created date
20240	001	Mahmood_Main	aaaaa		0	0	0	Yes	Yes	07-05-2021
20234	001	Child_lab_inova_M...	Lahore		0	0	0	Yes	Yes	07-05-2021
20229	JL	Jinah lab	Lahore		10	10	2000	No	Yes	19-04-2021
20197	Doc12	Doctors lab	Lahore		10	10	100	No	Yes	24-03-2021
20189	Hum	Hum CC	Lahore		10	10	1000	No	Yes	19-03-2021
20184	Chu00	Chugzal	Lahore		10	30	1000	No	Yes	17-03-2021
20177	jinah	Jinnah road	Lahore		10	10	0	No	Yes	15-03-2021
20168	001	Inova lab_Main	DHA		0	0	11000	Yes	Yes	07-05-2021

Showing 25 records per page. 1-8 Of 8

<<Start < 1 > End>>

Add Collection Centre

Collection centers are the outlets where a patient's samples are collected for tests. There can be laboratories that have multiple collection centers across a town but tests are performed in one central laboratory. To add a collection center, use the 'Collection center' menu and select the 'Add Collection Centre' option. Fill in the form with the collection center's details.

Set discount:

- The lab can set discounts per test and per day.
- Can set a daily counter.

ONELABz Dashboard Sample collection Test result Patient Accounting Stock Quality control Reports Settings Access control Admin Inova lab Inova lab

Add collection centre [← Back to collection centres listing](#)

Collection centre

Collection center code *
Code

Collection center name *
Name

Address *
Address

Maximum discount per test (%)
0

Maximum discount per day (%)
0

Daily counter start value
0

Head office

Enabled

Lab *
Inova lab

Description
Description

[Save](#)

Collection Centres

All collection centers which are added in OneLABz can be viewed on this screen.

ONELABz Dashboard Sample collection Test result Patient Accounting Stock Quality control Reports Settings Access control Admin Inova lab Inova lab

Lab settings [+ Add new collection center](#)

Collection center Employee Doctor

Export
Copy CSV Excel PDF Print

Search by <--Please select a value--> Search value Search value [Search](#) [Reset](#)

Collection centre id	Code	Collection center	Address	Description	Max disc per test (%)	Max disc per day (%)	Counter start	Head office	Enabled	Created date
20240	001	Mahmood_Main	aaaaa		0	0	0	Yes	Yes	07-05-2021
20234	001	Child_Lab_Inova_M...	Lahore		0	0	0	Yes	Yes	07-05-2021
20229	JL	Jineh lab	Lahore		10	10	2000	No	Yes	19-04-2021
20197	Doc12	Doctors lab	Lahore		10	10	100	No	Yes	24-03-2021
20189	Hum	Hum CC	Lahore		10	10	1000	No	Yes	19-03-2021
20184	Chu00	Chughtai	Lahore		10	30	1000	No	Yes	17-03-2021
20177	Jinah	Jinnah road	Lahore		10	10	0	No	Yes	15-03-2021
20168	001	Inova lab_Main	DHA		0	0	11000	Yes	Yes	07-05-2021

Showing 25 records per page. 1-8 Of 8

«Start < 1 > End»

Activate Windows
Activate Windows

Adding Employees

To add an employee, use the 'Employee' menu and select the 'Add Employee' option. Fill in the form with the employee's details.

The screenshot shows the 'Add employee' form in the OneLABz system. The form is divided into two main sections: 'Employee' and 'User details'. The 'Employee' section contains various fields for personal and professional information, including First name, Last name, Phone number, CNIC, Address, Gender, Age, Joining date, Salary (Rs.), Father/husband name, Email, Marital status, Lab, and Work collection center. A checkbox for 'Allow web login' is checked. A red box highlights the 'Marketing representative' checkbox. The 'User details' section includes fields for User name (Email), Password, and Employee role. A 'Save' button is located at the bottom left of the form. A watermark 'Activate Windows' is visible at the bottom right of the screenshot.

If user want to web access to the employee, click on the 'Allow web login' checkbox. Fill in the employee's login credentials by providing a 'Username' and 'Password'.

Click the 'Save' button to store employee details.

Note: there are three types of default system roles in OneLABz. Admin (users who can administer all aspects of software); Pathologist (users who approve test results, a test's report can not be printed unless it is approved by a pathologist); Lab operator (users who collect test samples and print results after a pathologist's approval)

The screenshot shows the 'User details' section of the 'Add employee' form. It includes the 'Allow web login' checkbox, which is checked. Below it are the 'Login id' and 'Password' fields, both with 'Email' and 'Password' entered respectively. A 'Save' button is located at the bottom right of this section. A watermark 'Copyright ©2020 Powered by TechXect' is visible at the bottom of the screenshot.

On the Employee screen, a checkbox is 'Marketing representative'.

A lab set the target on MR employees.

Collect sample referred doctors.

Employees Listing

Admin and Pathologist users can see the list of Employees. Employees can be searched using the 'Search by' criteria; to clear the search filter the 'Reset' button can be used. To copy records in the clipboard, export them in CSV, Excel, or PDF files or print them, users can choose from the appropriate export buttons displayed on the top left of the listing. The number of records shown per page can be altered from the paging dropdown at the botto left

of the page. Arrows on the bottom right of the page can be used to easily navigate through multiple pages. Records can be sorted by clicking column headers.

Pathologist

Pathologists approve test results. Admins or lab operators cannot print or send test results to doctors or patients without approval from the pathologist.

Lab settings

Collection center: **Employee** Doctor

Export: Copy CSV Excel PDF Print

Search by: <--Please select a value--> Search value: Search value Search Reset

Employee id	First name	Last name	Phone No.	Address	Collection center	Email	Created date	Created by
10414	Xcvcvxcv	Cvcvcvcx	324324324	dfgfdg	Mahmood_Main	ffdfgfdg@gmail.com	21-04-2021	MahmoodMahmood
10413	Mahmood	Mahmood	564654	aaaaa	Mahmood_Main	dsdsdsds@gmail.com	21-04-2021	Inova labInova lab
10394	Child_lab_inova	Child_lab_inova	65465654	Lahore	Child_lab_inova_M...	abc@gmail.com	20-04-2021	Inova labInova lab
10388	Employee	Testing	65165465	Lahore	Inova lab_Main	emply@onelabz.com	19-04-2021	Inova labInova lab
10336	Sania	Hussan	5464654	Lahore	Doctors lab	saa@onelabz.com	24-03-2021	Inova labInova lab
10333	Mr-ahsan	Ali	6554646	Lahore	Hum CC	ahs@onelabz.com	19-03-2021	Inova labInova lab
10330	Mr	Huanan	45465	Lahore	Chugtai	hunain@onelabz.com	18-03-2021	Inova labInova lab
10329	Mr	Hussan	654654654	Lahore	Chugtai	hus@onelabz.com	18-03-2021	Inova labInova lab
10326	CC	3	3654	Lahore	Jinnah road	CC3@onelabz.com	16-03-2021	Inova labInova lab
10320	CC	2	654512	Lahore	Jinnah road	CC2@onelabz.com	15-03-2021	Inova labInova lab
10319	Jinnah	CC	654654	Lahore	Jinnah road	laa@onelabz.com	15-03-2021	Inova labInova lab
10317	Inova	Receptionist	54654	Lahore	Inova lab_Main	inova_Rec@onelabz.com	15-03-2021	Inova labInova lab
10297	Inova	Pathologist	042+85+652+6	lahore	Inova lab_Main	path@onelabz.com	03-03-2021	Inova labInova lab
10296	Inova	Lab operator	04265566	Lahore	Inova lab_Main	la@onelabz.com	03-03-2021	Inova labInova lab
10295	Inova lab	Inova lab	58987654	DHA	Inova lab_Main	farzana.akram@xtsol.com	03-03-2021	SuperAdmin

Showing 25 records per page. 1-15 Of 15

Performed by: Muneer Hassan

Performed on date: 30/03/2020 17:10

Overall result: Overall result

Tests

Parameter	Result	Result in words	Abnormal	Male	Female	Child	Child2	Infant	New born	Min max value in words	Unit	Parameter type
Blood	0	abnorm	<input checked="" type="checkbox"/>	1-2	1-2	1-2	1-2	1-2	1-2	normal	mm/l/g	Che

Test is approved by: Ali Ahmad

Notes

Comments

Method

Print result Resend email to patient Resend email to doctor Multiple prints Multiple prints on letter head

Lab Operator

Lab operators have fewer access rights as compared to admin or pathologists, and they cannot print a test result unless a pathologist approves it.

Adding Doctors

OneLABz stores doctors' details in the system to keep track of doctors who refer patients for test results, and to calculate their referral share percentage. To add a doctor, use the 'Doctor' menu and select the 'Add Doctor' option. Fill in the form with the doctor's details.

Add doctor

Lab: Inova lab | Collection center: Inova lab_Main | Referred by employee: <--Please select a value-->

First name: First Name | Last name: Last Name | Clinic name: Clinic Name

Phone number: Phone number | CNIC: CNIC | Address: Address

Regular share percentage: 0 | Gender: <--Please select a value--> | Email: Email

Qualification: Qualification | Print on report | Print order (Your limit for doctor on report is:5. Remaining limit is: -9) 0

Always send Email | Always send SMS | Always send WhatsApp

Allow web login

User details: User name: Email | Password: Password | Doctor role: Inova Lab_Doctor

[Back to doctors listing](#)

User can set the Number of doctors which user wants to print on report.

Doctors

All doctors who are added in OneLABz can be viewed on this screen.

Lab settings

Collection center | Employee | **Doctor**

Export: Copy | CSV | Excel | PDF | Print

Search by: <--Please select a value--> | Search value: Search value | Search | Reset

Doctor id	Doctor number	First name	Last name	Phone No.	Address	Email	Created date	Created by	Modified
20250	DR20210414-0-001-00...	Asdsd	Dstdsf	865231	0		14-04-2021	Inova lab Inova lab	14-04-2021
20214	DR20210324-1-Doc12...	Dr	Sana	6544984	Lahore		24-03-2021	Inova lab Inova lab	14-04-2021
20213	DR20210324-1-Hum-0...	Dr	Balkees	6565464	Lahore		24-03-2021	Inova lab Inova lab	24-03-2021
20212	DR20210324-1-Hum-0...	DR	Zohaib	6566556	Lahore		24-03-2021	Inova lab Inova lab	24-03-2021
20211	DR20210324-1-Hum-0...	DR	Saaim	654484	Kasur		24-03-2021	Inova lab Inova lab	24-03-2021
20205	DR20210319-1-Hum-0...	Dr	Israr ul haq	65646544	Lahore		19-03-2021	Inova lab Inova lab	19-03-2021
20204	DR20210319-1-Hum-0...	Dr	Salif	5156554	Lahore		19-03-2021	Inova lab Inova lab	19-03-2021
20202	DR20210318-1-Chu00...	Dr	Shakil	2465646	Lahore	farzana.akram@txtsol.com	18-03-2021	Inova lab Inova lab	28-04-2021
20200	DR20210316-1-001-00...	Dr	Sara	654	Lahore		16-03-2021	Inova lab Inova lab	16-03-2021
20198	DR20210316-1-jinah-0...	DR	Tahira	456	Lahore		16-03-2021	Inova lab Inova lab	16-03-2021
20197	DR20210316-1-jinah-0...	Dr	Pankraj	5468	Lahore		16-03-2021	Inova lab Inova lab	16-03-2021
20196	DR20210316-1-jinah-0...	Dr	Rashid	4562	Lahore		16-03-2021	Inova lab Inova lab	16-03-2021
20194	DR20210315-1-001-00...	Dr	Shahryar	064565	Lahore		15-03-2021	Inova lab Inova lab	15-03-2021
20192	DR20210312-1-001-00...	Dr	Shakoore	04533593	Lahore		12-03-2021	Inova lab Inova lab	12-03-2021
20188	DR20210312-1-001-00...	DR	Nasir	354654	Lahore		12-03-2021	Inova lab Inova lab	12-03-2021
20187	DR20210312-1-001-00...	Dr	Shazia	01246	Lahore		12-03-2021	Inova lab Inova lab	12-03-2021
20185	DR20210312-1-001-00...	Dr	Shahina	65654654654	Lahore		12-03-2021	Inova lab Inova lab	12-03-2021

Add Lab

To add a lab, use the 'Admin' menu and select the 'Add Lab' submenu. Fill in the form with the 'Lab', 'Login details', and assign Lab features.

Add lab details:

Dashboard Sample collection Test result Patient Expense Stock Quality control Reports Settings Access control Admin

Add lab

Lab

Lab name *
Lab Name

Parent lab *
Lab for user manual

Slogan *
Slogan

Address line 1 *
Address line 1

Address line 2
Address line 2

Address line 3
Address line 3

Landline number *
Landline number

Mobile number *
Mobile number

WhatsApp number *
WhatsApp number

Website *
Website

Email *
Email

Date format *
<--Please select a value-->

Currency *
<--Please select a value-->

Activate Wi
Go to Settings

Add Login details:

Allow web login

Login details

Login id *
Login id

Password *
Password

Save

Assign Features:

Lab features

Available features

- Menus
- Dashboard
- Sample collection
- Collect sample
- Samples
- Test result
- Add test result
- Test results
- Patient
- Add patient

Selected features


<< < > >>

Save

Labs


Admin and Pathologist users can see the list of Labs. Employees can be searched using the 'Search by' criteria; to clear the search filter the 'Reset' button can be used. To copy records in the clipboard, export them in CSV, Excel, or PDF files or print them, users can choose from the appropriate export buttons displayed on the top left of the listing. The number of records shown per page can be altered from the paging dropdown at the bottom left of the

page. Arrows on the bottom right of the page can be used to easily navigate through multiple pages. Records can be sorted by clicking column headers.



Home sample collect request

Without login book free home sample by phone call.



Add home sample collection request

By using the add home sample collection request screen can add samples, fill all mandatory fields and change status.

Home sample collection request

Lab *
Sahiwal Lab

Collection center *
Sahiwal Lab_Main

Patient

First name *

Last name *

Phone number *

Family number

CNIC

Address *

Patient type *
--Please select a value--

Gender
--Please select a value--

Blood group
--Please select a value--

Date of birth

Age in day(s) - month(s) - year(s) *

Father/husband name *

Email *

Marital status *
--Please select a value--

Passport number

Allow web login

Home sample collection detail

Collection request status *
New

Collection request source
--Please select a value--

Preferred collection time *

Collection address *

Scheduled collection time

Tests to perform *

Who will collect

Home Sample

All home samples which are added in OneLABz can be viewed on this screen.

Home sample collection requests

Export
Copy CSV Excel PDF Print

Search by --Please select a value-- Search value Search Reset


Id	Patient	Phone no.	Email	Address	Tests to perform	Status	Source	T/R No.
106	Patient 2	6546554654	patinet2@onelabz...	Lahore	DNA,MRI	Complete	Phone call	TR20210430-0-001-0000
105	Patient 1	654655	emall@email.com	Lahore	Blood test	Collection scheduled	Online	
92	Khalid Ayub	654654	344@gmail.com	Lahore	blood	Complete	Online	TR20210421-0-001-0016
91	Irfan Ahmad	0302562289	syeryr@email.com	Lahore	CBC	Sample collected	Online	
90	Ayza Khan	0302562289	ayza@email.com	Lahore	Blood	Complete	Online	TR20210421-0-001-0015
89	Iqra Shokat	65465564	saira@OneLABz.com	Multan	CBC	Complete	Online	TR20210421-0-001-0005
88	Sam Ben	654654654654	marcus@OneLABz...	England	CBC	Complete	Online	
87	Malik Shahzain	6565654	shahzain@OneLAB...	Multan	CBC, Bacterial infection	Complete	Online	TR20210421-0-001-0003
85	Shahid Iqbal	65454654	shahid@OneLabz...	Lahore	CBC,Blood test	Sample collected	Online	
84	Patient_3_sahiwal_Lab	0302562289	p22@email.com	Lahore	CBC,MRI	Sample collected	Online	
83	Irfan Sahiwal lab	0302562289	patient2@email.com	Lahore	CBC,Blood	Sample collected	Online	

Showing 25 records per page. 1-11 Of 11

«Start < 1 > End»

Custom notification


If you want to send email or SMS notification use this screen to fill all mandatory fields.



Admin settings

Add Currency

To 'add currency', use the 'Admin Settings' menu and select the 'Add currency' tab. Fill in the form with the currency details.



Currencies

Admin users can see the list of Currencies. Currencies can be searched using the 'Search by' criteria; to clear the search filter the 'Reset' button can be used. To copy records in the clipboard, export them in CSV, Excel, or PDF files or print them, users can choose from the appropriate export buttons displayed on the top left of the listing. The number of records

shown per page can be altered from the paging dropdown at the bottom left of the page. Arrows on the bottom right of the page can be used to easily navigate through multiple pages. Records can be sorted by clicking column headers.

Admin settings

Measuring unit

Export: Copy, CSV, Excel, PDF, Print

Search by: <--Please select a value--> Search value: Search value Search Reset

Currency id	Name	Code	Symbol	Deleted	Created date	Lab
13	Us_dollar_d	12	\$\$	No	02-03-2021	New lab for sp testing
12	Currency	Cur	#	Yes	22-02-2021	New lab for sp testing
11	PKR	001	Rprt	No	16-02-2021	New lab for sp testing
10	Tm	12s	@	Yes	29-12-2020	Regression
9	Fake Currency	Fak	Fake Curre	Yes	28-12-2020	After major change
8	MNO	01	PKR 123	Yes	18-11-2020	Admin
7	JKL	21	PKR	Yes	20-10-2020	Lab for user manual
6	GHI	010	RS	Yes	01-07-2020	Admin
5	DEF	121	RS.2	Yes	29-06-2020	Shifa Lab
4	ABC	02	RS.7	Yes	25-06-2020	Hormone Lab
3	Euro	222	€	No	12-05-2020	Lab2
2	Dollar	123	\$	No	12-05-2020	Lab2
1	Rupee	101	Rs.	No	12-05-2020	Lab2

Showing 25 records per page. 1-13 Of 13

Add measuring unit

To add a measuring unit, use the 'Admin' menu and select the 'Add measuring unit' tab. Fill in the form with the measuring unit's details.

Add measuring unit

Measuring unit

← Back to measuring units listing

Name *
Name

Measuring type *
<--Please select a value-->

Abbreviation *
Abbreviation

Smallest unit

Description
Description

Save

Measuring units

Admin users can see the list of Measuring units. Measuring units can be searched using the 'Search by' criteria; to clear the search filter the 'Reset' button can be used. To copy records in the clipboard, export them in CSV, Excel or PDF files or print them, users can choose from the appropriate export buttons displayed on the top left of the listing. The number of records shown per page can be altered from the paging dropdown at the bottom left of the page. Arrows on the bottom right of the page can be used to easily navigate through multiple pages. Records can be sorted by clicking column headers.

ONE LABZ Dashboard Sample collection Test result Patient Accounting Stock Quality control Reports Settings Access control Admin Inova lab Inova lab

Admin settings + Add new measuring unit

Currency Measuring unit Items

Export

Search by <--Please select a value--> Search value Search value

Measuring unit id	Name	Measuring type	Abbreviation	Smallest unit	Description	Created date	Lab
17	Gass	Gass	gess	Yes		25-03-2021	Computer
16	Liter	Liquid	l	Yes	liter	04-03-2021	Computer
15	Kilo	Mass	k	Yes	kg	04-03-2021	Computer
14	Garm	Mass	G1	Yes		04-01-2021	Admin
13	H2O	Plasma	water	Yes		04-01-2021	Lab_Test_Packiges
12	Solida	Mass	s	Yes	aaa	29-12-2020	Regression
11	Envelop	Mass	Envelop	Yes		18-12-2020	Lab_Test_Packiges
10	Test	Mass	sr	Yes		18-11-2020	Admin
9	Syringe	Mass	sr	Yes	Syringe	20-10-2020	Lab for user manual
8	MG	Mass	miligram	No		29-06-2020	Shifa Lab
7	Kilo gram	Liquid	xyz	Yes		25-06-2020	Hormone Lab
6	Ml	Liquid	hgfoib	Yes	jhgvihjv	22-06-2020	Excel Lab
5	Piece	Mass	Piece	Yes		22-06-2020	Excel Lab
4	Litter	Mass	Piece	Yes	Piece	12-05-2020	Lab2
3	g/dl	Gass	g/dl	Yes	Gram per decilitre	12-05-2020	Lab2
2	Abc	Liquid	g	Yes	Gram	12-05-2020	Lab2
1	Milli Gram	Mass	mg	Yes	Milli Gram		

Activate Windows
 Go to Settings to activate Windows.

Add item

To record inventory and availability of chemicals that are used for different tests, the chemicals are added to the system. The inventory of chemicals is automatically reduced when a test is performed in which the chemical is used. To add a chemical, use the 'Admin chemical' menu and select the 'Add admin chemical' option. Fill in the form with the chemical's details and use the 'Save' button to store the input values.

ONE LABZ Dashboard Sample collection Test result Patient Accounting Stock Quality control Reports Settings Access control Admin Inova lab Inova lab

Add item ← Back to items listing

Item

Name * Description

Stock category * <--Please select a value-->

Alert threshold * Active

Threshold measuring type * <--Please select a value-->

Measuring type * <--Please select a value-->

- <--Please select a value-->
- Mass
- Liquid
- Gass
- Plasma

Items

All chemicals which are added in OneLABz can be viewed on this screen.

ONE LABS Dashboard Sample collection Test result Patient Accounting Stock Quality control Reports Settings Access control Admin Inova lab Inova lab

Admin settings + Add new item

Currency Measuring unit **Items**

Export Copy CSV Excel PDF Print

Search by <-Please select a value--> Search value Search value Search Reset

Item id	Name	Description	Measuring type	Stock category	Threshold	Expiry day(s) threshold	Created date	Created by	Modified date	Mo
10476	Chlorine Liquid	Chlorine Liquid	Liquid	Refrigerator Chemical	10	2	28-04-2021	Inova lab Inova lab	28-04-2021	Inov
10471	Phosphors	Phosphors	Mass	Refrigerator Chemical	0	0	08-04-2021	Inova lab Inova lab	07-05-2021	
10470	HCL	HCL	Liquid	Refrigerator Chemical	0	0	24-03-2021	Inova lab Inova lab	07-05-2021	
10469	Nitrogen	Nitrogen	Gass	Other Items	0	0	19-03-2021	Inova lab Inova lab	07-05-2021	
10468	Florine	Florine	Liquid	Refrigerator Chemical	0	0	19-03-2021	Inova lab Inova lab	07-05-2021	
10467	Chlorine	Chlorine	Liquid	Refrigerator Chemical	0	0	19-03-2021	Inova lab Inova lab	07-05-2021	
10466	Hydrogen	Hydrogen	Mass	Devices / Strips	0	0	19-03-2021	Inova lab Inova lab	07-05-2021	
10465	Yellow Tips	Yellow Tips	Mass	Devices / Strips	0	0	15-03-2021	Inova lab Inova lab	07-05-2021	
10464	Cover slip	Cover Slip	Mass	Devices / Strips	0	0	03-03-2021	Inova Lab operator	07-05-2021	
10463	Blood Testing kit	xyz	Mass	Other Items	5	15	01-07-2020	Super Admin	07-05-2021	
10462	Glass tube	Glass Tube	Mass	Devices / Strips	0	0	17-03-2021	Inova lab Inova lab	07-05-2021	
10452	Phosphors	Phosphors	Mass	Refrigerator Chemical	0	0	08-04-2021	Inova lab Inova lab	07-05-2021	
10451	HCL	HCL	Liquid	Refrigerator Chemical	0	0	24-03-2021	Inova lab Inova lab	07-05-2021	
10450	Nitrogen	Nitrogen	Gass	Other Items	0	0	19-03-2021	Inova lab Inova lab	07-05-2021	
10449	Florine	Florine	Liquid	Refrigerator Chemical	0	0	19-03-2021	Inova lab Inova lab	07-05-2021	
10448	Chlorine	Chlorine	Liquid	Refrigerator Chemical	0	0	19-03-2021	Inova lab Inova lab	07-05-2021	
10447	Hydrogen	Hydrogen	Mass	Devices / Strips	0	0	19-03-2021	Inova lab Inova lab	07-05-2021	
10446	Yellow Tips	Yellow Tips	Mass	Devices / Strips	0	0	15-03-2021	Inova lab Inova lab	07-05-2021	
10445	Cover slip	Cover Slip	Mass	Devices / Strips	0	0	03-03-2021	Inova Lab operator	07-05-2021	

Notification settings

For notification use these screens

- Email notification
- SMS notification

Email notification

For more addition click on the New add button

ONE LABS Dashboard Sample collection Test result Patient Accounting Stock Quality control Reports Settings Access control Admin Sahiwal Lab Sahiwal Lab

Notification settings + New email template

Email templates Message templates

Export Copy CSV Excel PDF Print

Search by <-Please select a value--> Search value Search value Search Reset

Id	Email template	Subject	Created date	Created by	Modified date	Modified by	Lab
1942	ForgotPassword	Instructions to reset OneLABs password	2021-02-03	Super Admin	2021-02-03	Super Admin	Sahiwal Lab
1941	CreateLabNotification	Your lab is created	2020-03-27	Super Admin	2020-03-27	Super Admin	Sahiwal Lab
1940	SendChemicalExpiryNotification	SendChemicalExpiryNotification	2020-02-25	Super Admin	2020-02-25	Super Admin	Sahiwal Lab
1939	ChemicalThresholdNotification	Chemical Stock Threshold Notification	2020-01-04	Super Admin	2020-01-04	Super Admin	Sahiwal Lab
1938	CreateDoctorNotification	Welcome To Rahbar Labs Family	2020-01-04	Super Admin	2020-01-04	Super Admin	Sahiwal Lab
1937	PatientSampleCollectedModifyNotification	Sample Collection Update	2019-10-24	Super Admin	2019-10-24	Super Admin	Sahiwal Lab
1936	PatientSampleCollectedWithCompletePay...	Thank You For Visiting Rahbar Labs	2019-10-24	Super Admin	2019-10-24	Super Admin	Sahiwal Lab
1935	CreatePatientNotification	Welcome To Rahbar Labs Family	2019-10-24	Super Admin	2019-10-24	Super Admin	Sahiwal Lab
1934	DoctorNotifyPatientTestReport	Patient Test Report Rahbar Labs	2019-10-22	Super Admin	2019-10-22	Super Admin	Sahiwal Lab
1933	DoctorNotifySampleCollected	Patient Sample Collected At Rahbar Labs	2019-10-22	Super Admin	2019-10-22	Super Admin	Sahiwal Lab
1932	PatientSampleCollectedWithDuePayment	Thank You For Visiting Rahbar Labs	2019-07-08	Super Admin	2019-07-08	Super Admin	Sahiwal Lab
1931	PatientNotifyTestReportWithCompletePa...	Test Report	2019-07-07	Super Admin	2019-07-07	Super Admin	Sahiwal Lab
1930	PatientNotifyTestReportWithDuePayment	Pay Your Pending Dues	2019-07-17	Super Admin	2019-07-17	Super Admin	Sahiwal Lab

Showing 25 records per page. 1-13 Of 13

SMS notification

For more addition click on the New add button

Notification settings

Export: Copy, CSV, Excel, PDF, Print

Search by: <--Please select a value--> Search value: Search value Search Reset

Id	Message template	Created date	Created by	Modified date	Modified by	Lab
No record found.						

Showing 25 records per page. 0-0 Of 0

Settings

Collection center rate list:

The Collection center sets rate against any test, commission, discount, and copy that rate into other collection centers.

Collection centre rate list

Lab: Sahiwal Lab Collection centre: Sahiwal_Collection_Center Test category: Bacteriology

Discount %: 0 Commission (Rs.): 0

Apply Copy rate list

Test name	Lab rate (Rs.)	Collection centre rate (Rs.)	Commission (Rs.)	Discount bear by collection centre(%)	Discount bear by lab(%)	
Bacter	500.000	0.000	0.000	0.000	100.000	Update

Save

Lab Configurations

Admin users can see the list of Lab configurations. Employees can be searched using the 'Search by' criteria; to clear the search filter the 'Reset' button can be used. To copy records in the clipboard, export them in CSV, Excel, or PDF files or print them, users can choose from the appropriate export buttons displayed on the top left of the listing. The number of records shown per page can be altered from the paging dropdown at the bottom left of the page. Arrows on the bottom right of the page can be used to easily navigate through multiple


Lab configurations

Export: Copy, CSV, Excel, PDF, Print

Search by: Lab configurations Search value: Search value Search Reset

Configuration id	Lab name	Configuration name	Description	Value
434	Madina lab	Patient due summary report path	PatientDueSummaryReport	../PrintTemplate/MadinaLab/PatientDu
433	Madina lab	Multiple test result report without header p...	MultipleTestResultReportWithoutHeaderPath	../PrintTemplate/MadinaLab/MultipleTe
432	Madina lab	Multiple test result report path	MultipleTestResultReportPath	../PrintTemplate/MadinaLab/MultipleTe
431	Madina lab	Single test result report path	SingleTestResultReportPath	../PrintTemplate/MadinaLab/TestResult
430	Madina lab	Summary reports path	SummaryReportsPath	../PrintTemplate/GeneralReportPatienti


pages. Records can be sorted by clicking column headers. Admin can modify configuration values



Access control

Add roles

When adding any lab by default it creates five different types of roles, 'Admin', 'Pathologist', Lab operator, Doctor and patient, if admin assign features but for 'Pathologist', Lab operator, Doctor and patient have limited features.



Admin has different types of role:

Admin

Admin has all roles in which all graphs, access rights.

Lab operator:

Lab operators have fewer access rights as compared to Admin.

Pathologist:

Pathologists have just approval rights.

Doctor:

Doctor have access rights

- Collect sample
- Samples
- Test results
- Doctor share report

Patient:

Patients just can view their test results.

- Test result listing.

Collection center:

Collection center have access rights

- Collect sample
- Samples

Reception:

Have access rights

- Collect sample
- Samples listing

Marketing Representative:

Market representative have access rights

- Dashboard
- Lab settings
- Doctor
- Add doctor

Roles

Admin users can see the list of roles. Roles can be searched using the 'Search by' criteria; to clear the search filter the 'Reset' button can be used. To copy records in the clipboard, export them in CSV, Excel, or PDF files or print them, users can choose from the appropriate export buttons displayed on the top left of the listing. The number of records shown per page can be altered from the paging dropdown at the bottom left of the page. Arrows on the

bottom right of the page can be used to easily navigate through multiple pages. Records can

Role id	Name	Lab name	Description	Active	Created date
137	Madina lab_Patient	Madina lab	Madina lab_Patient	Yes	20/10/2020
136	Madina lab_Doctor	Madina lab	Madina lab_Doctor	Yes	20/10/2020
135	Madina lab_LabOperator	Madina lab	Madina lab_LabOperator	Yes	20/10/2020
134	Madina lab_Pathologist	Madina lab	Madina lab_Pathologist	Yes	20/10/2020
133	Madina lab_Admin	Madina lab	Madina lab_Admin	Yes	20/10/2020

be sorted by clicking column headers.

Patient

Sign-inAdd patient

To add a patient, use the 'Patient' menu and select the 'Add Patient' option. Fill in the form with the patient's details.

First name *

Last name *

Phone number *

CNIC

Address *

Patient type *

Age in day(s) - month(s) - year(s) *

Father/husband name *

Email

Marital status *

Check the 'Allow web login' checkbox, if web access is allowed for a patient.

Allow web login

User details

Login id *

Email

Password *

Password

Save

Activate Windows

Go to Settings to activate Windows.

Copyright ©2020 Powered by **TechXect**

Patients Listing

All patients who are added in OneLABz can be viewed on this screen.

Patient id	Patient first name	Patient last name	Patient phone#	Patient address	Patient email	Created date
40488	Kaleem	Alli	0343406...	zarafarm pakpattan chok sahiwal.	sarmad.abbas@txtsol.com	1/12/2020
40550	Hajira	Bibi	0300416...	wapda sahiwal.	sarmad.abbas@txtsol.com	1/15/2020
40602	Shoukat	Mumtaz	0300792...	wapda sahiwal.	sarmad.abbas@txtsol.com	1/16/2020
50839	Mohammad	Anees	0306693...	wapda colony	sarmad.abbas@txtsol.com	1/22/2020
40600	Javeria		0300969...	wahdat colony swl	sarmad.abbas@txtsol.com	1/16/2020
40309	Abdul	Waheed	0301693...	wahdat colony arif road Sahiwal.	sarmad.abbas@txtsol.com	1/9/2020

Stock

Manage Stock

To load or dispose of chemicals' stock, the 'Manage stock' screen is used. The inventory of chemicals is automatically reduced when a test is performed, in which the chemical is used. To add or dispose of a chemical, use the 'Stock' menu and select the 'Manage stock' option. Fill in the form with the stock details and use the 'Save' button to store the input values.

Manage stock

Batch number *

Chemical *

Measuring unit *

Quantity *

Stock action *

Issue date *

Expiry date

Buying date

Comment

Save

Stock Listing

All stocks which are added in OneLABz can be viewed on this screen.

Stock

Export

Copy CSV Excel PDF Print

Search by <--Please select a value--> Search value Search value Search Reset

Stock id	Chemical	Measuring unit	Available quantity	Used quantity	Total quantity	Buying date	Issue date
10075	Washing Brush	Piece	10		10		
75	Acetic Acid	MilliLiter	297.85		300		
74	Envelop	Piece	67		1400		
73	Letter Pad	Piece	963		2300		
72	EDTA Liquid	MilliLiter	856.8		900		
71	lyse	MilliLiter	406		950		
70	Diluent	MilliLiter	6401		26423		
69	Sani Plast	Piece	523		1764		
68	Surgical Gloves	Piece	80	0	80		
67	Syringe Stoper	Piece	38	962	1000		
66	Serum Cups	Piece	306	14	320		

Activate Windows
Go to Settings to activate Windows.

Quality Control

Add Quality Control values

Input all mandatory fields like (Tests, methods, Measuring Unit, Min value, Target value, Max value, created date) and click the 'save' button to save.

Dashboard Sample collection Test result Patient Expense Stock Quality control Reports Settings Access control Admin

admin admin

Add quality control values

Add quality control values
Quality controls listing

Quality control values

Test * <--Please select a value--> Method * <--Please select a value--> Measuring unit * Measuring unit

Min value * 0 Max value * 0 Target value * 0

Description Description Start date * 23/10/2020 10:47

Save

Quality control listing

All quality control system records which are added in OneLABz can be viewed on this screen.

Dashboard Sample collection Test result Patient Expense Stock Quality control Reports Settings Access control Admin admin admin

Quality controls listing

Add quality control values
Quality controls listing

Export
Copy CSV Excel PDF Print

Search by <--Please select a value--> Search value Search value Search Reset

Quality control id	Test	Method	Description	Start date	Inspection date
13	Albumin	DGKC IFCC 1994 37 °C		20/04/2020	
12	Albumin	DGKC IFCC 1994 37 °C		20/04/2020	
11	Total Protein	Biuret		18/01/2020	
10	Uric Acid	Uricase		18/01/2020	
9	Urea	Urease UV		18/01/2020	
8	Creatinine	Jaffe Without Compensation		18/01/2020	
7	Cholesterol	CHOD-PAP		18/01/2020	
6	Calcium Level	Phosphonazo		18/01/2020	
5	SGOT (AST)	IFCC With P5P 37 °C		18/01/2020	
4	SGPT(ALT)	IFCC With P5P 37 °C		18/01/2020	
3	Alkaline Phosphatase (IFCC 37 °C		18/01/2020	

Sample Collection

Collect Sample

In order to record tests that should be conducted for a patient, use the 'Add sample collection' screen. Use the 'Sample collection' menu and select 'Add sample collection'.

The patient for whom the test is performed should be selected in the 'Patient' field. Use the 'Add patient' button if a patient does not already exist in the system.

Appropriate checkboxes for email, SMS, or WhatsApp can be selected if a patient or doctor needs to be notified when test results are ready.

Dashboard Sample collection Test result Patient Expense Stock Quality control Reports Settings Access control Admin admin admin

Collect sample
Samples

Other detail

External test reference number
External test reference number

Referred by doctor
Doctor

Doctor clinic
Doctor clinic

Flight details

Airline *
Airline

Flight number
Flight number

Flight date *
2021-05-17 12:40

Ticket number
Ticket number

Origin city
Origin city

Destination city
Destination city

PNR number
PNR number

Comments
Comments

Address
Address

Marital status
<--Please select a value-->

Collection center *
Main lab

Collection date *
23/10/2020 10:51

Referred by doctor
Doctor

Comments
Comments

If a patient takes flight that detail is also given on the collect sample page. Check the checkbox and fill in the flight details.

Tests to perform

Tests that are to be performed for the patient should be added in the 'Test' section by using the 'Add test' button.

Code	Category	Name	Specimen	Specimen collection type	Charges (Rs.)	
<input type="text" value="Code"/>	<input type="text" value="Select a value"/>	<input type="text"/>	<input type="button" value="Remove"/>			

Use the 'Remove' button at the end of the record to delete a particular test record.

Sample status

Admin can set the sample status, if the test report requires urgent to check the checkbox.

Then its status became urgent.

Code	Category	Name	Specimen	Specimen collection type	Sample status	Urgent	Charges (Rs.)
<input type="text" value="Code"/>	<input type="text" value="Category"/>	<input type="text" value="Name"/>	<input type="text" value="Specimen"/>	<input type="text" value="Select a value"/>	<input type="text" value="Select a value"/>	<input type="checkbox"/>	<input type="text"/>

Notification

Checked checkboxes for notification. if patient or doctor wants emails then check email checkboxes for SMS and WhatsApp notification select SMS or WhatsApp checkboxes

Notification

Email to patient SMS to patient WhatsApp to patient

Email to doctor SMS to doctor WhatsApp to doctor

Payment details

Fill in the payment details in the 'Payment details section. Please note that the standard configured charges are applied by default for the tests.

To apply a discount use the 'Discount (%)' or 'Discount amount' fields.

If you wish to print an invoice receipt, use the 'Save and print invoice' button. Use the 'Save' button to just store the details in the system.

Payment details

Gross amount (Rs.) *	Discount (%) - [Allowed discount 0%]	Discount amount (Rs.)	Net amount (Rs.)
<input type="text" value="0.000"/>	<input type="text" value="0.00"/>	<input type="text" value="0.000"/>	<input type="text" value="0.000"/>
Paid amount (Rs.)	Doctor share (Rs.)	Due amount (Rs.)	Now paying (Rs.)
<input type="text" value="0.000"/>	<input type="text" value="0.000"/>	<input type="text" value="0.000"/>	<input type="text" value="Now paying"/>

Invoice print

After payment, an invoice is generated. On that print different type. Admin can print A4, A5, or thermal.

A4:

Modify collected sample 192.168.4.119

192.168.4.119/OneLabzWeb/TestInstanceGroup/PrintInvoice.html?It=3b8jn7TjO1no0DRvyryw29R8%2FY%2FT13NLIKVE1JzV3GKzmzRLN4aR%...
 Log in Webmail Login T-Mobile M2M Hub EdgeOS - ubnt Log in PTP - Admin QA - Google Drive ONE Deployment Q... Reading list

A4
 A5
 Thermal

ONE LABz

Print Print lab copy

Patient #: 40897
 Lab No.: 63343

TRACK ONLINE
 QR Code

Patient Name: Hamza Ali
 Father/Husband Name: Father
 Age/Sex: 21 (Y) - 8 (M) - 13 (D) / Male
 Marital status: Single
 Phone number: 923440041445
 Address: Cannt
 Passport No.:

M/R No.:
 Registered at:
 Sample collected on:
 Consultant:
 Email:
 CNIC:

CS-43
 PAF_Child_Lab_Main
 25-Jun-2021 12:28 PM
 PAF_Child_Lab PAF_Child_Lab
 Self
 hamza@email.com

Patient bill

S. No.	Test Name	Reporting Date Time	Rate
1	Bacteriology	26-Jun-2021 12:28 PM	3,000.00
2	blood test	25-Jun-2021 12:28 PM	300.00

VIEW LAB REPORTS ONLINE :
 To view report online, scan the above QR code.
 4 OR
 This report can be verified at OneLabz login, click on "View test result"

Test reference number 63343
 Patient Id 40897

Electronically verified receipt, no signature(s) required.

www.PAF.com , Tel: 65654 , Mobile: 65465 , 65654 , Cannt
 Powered by TechXect

Activate Windows
 Go to Settings to activate Windows.

A5:

Not secure | 192.168.4.119/OneLabzWeb/TestInstanceGroup/PrintInvoice.html?It=3b8jn7TjO1no0DRvyryw29R8%2FY%2FT13NLIKVE1JzV3GKzmzRLN4aR%...
 PTP - Admin OneLogin Log in Webmail Login T-Mobile M2M Hub EdgeOS - ubnt Log in PTP - Admin QA - Google Drive ONE Deployment Q... Reading list

A4
 A5
 Thermal

ONE LABz

Print Print lab copy

Patient #: 40897
 Lab No.: 63343

TRACK ONLINE
 QR Code

Patient Name: Hamza Ali
 Father/Husband Name: Father
 Age/Sex: 21 (Y) - 8 (M) - 13 (D) / Male
 Marital status: Single
 Phone number: 923440041445
 Address: Cannt
 Passport No.:

M/R No.:
 Registered at:
 Sample collected on:
 Consultant:
 Email:
 CNIC:

PAF_Child_Lab_Main
 25-Jun-2021 12:28 PM
 PAF_Child_Lab PAF_Child_Lab
 Self
 hamza@email.com

Patient bill

S. No.	Test Name	Reporting Date Time	Rate
1	Bacteriology	26-Jun-2021 12:28 PM	3,000.00
2	blood test	25-Jun-2021 12:28 PM	300.00
3	CBC	25-Jun-2021 01:28 PM	500.00
4	CH patho	29-Jun-2021 12:28 PM	2,500.00
5	DNA	05-Jul-2021 12:28 PM	2,300.00
6	Endocrinology	25-Jun-2021 01:28 PM	9,600.00
7	Hematology	26-Jun-2021 12:29 PM	2,300.00
8	Histo	26-Jun-2021 12:29 PM	2,500.00

VIEW LAB REPORTS ONLINE :
 To view report online, scan the above QR code.
 OR
 This report can be verified at OneLabz login, click on "View test result"

Test reference number 63343
 Patient Id 40897

Electronically verified receipt, no signature(s) required.

www.PAF.com , Tel: 65654 , Mobile: 65465 , 65654 , Cannt
 Powered by TechXect

Activate Windows
 Go to Settings to activate Windows.

23,000.00
 130
 300.00
 22,700.00
 23,000.00

Thermal invoice report:

● Thermal



Patient #: 
40897

Lab #: 
63343

PAF_Child_Lab

Cannt

Ph : 65654 , 65465

www.PAF.com

Name: Hamza Ali
 Gender: 21 (Y) - 8 (M) - 13 (D) / Male
 Date: 25-Jun-2021 12:28 PM
 M.R.#: 63343
 F/H Name: Father
 NIC:
 Address: Cannt
 Phone No: 923440041445
 Reference: PAF_Child_Lab_Main
 Consultant: Respected Dr.

Patient Bill

S#	Test Name	Rate
1	Bacteriology	3,000.00
2	blood test	300.00
3	CBC	500.00
4	CH patho	2,500.00
5	DNA	2,300.00
6	Endocrinology	9,600.00
7	Hematology	2,300.00
8

Samples

All sample collections which are added in OneLABz can be viewed on this screen.

The screenshot shows the 'Samples' page in the OneLABz application. At the top, there is a navigation bar with 'Dashboard', 'Sample collection', 'Test result', 'Patient', 'Accounting', 'Stock', 'Quality control', 'Reports', 'Settings', 'Access control', and 'Admin'. The 'Samples' page title is visible, along with a 'Samples settings' link. Below the title, there are 'Export' options (Copy, CSV, Excel, PDF, Print) and a search bar. The main content is a table with the following columns: T/R Id, T/R No., M/R No., Collection centre, Collection date, Patient, Patient phone#, Patient type / Age, and Doctor. The table contains 20 rows of sample data.

T/R Id	T/R No.	M/R No.	Collection centre	Collection date	Patient	Patient phone#	Patient type / Age	Dr
62552	TR20210517-0-001-0000		Sahiwal Lab_Main	2021-05-17	Muneeba Kafayat	9254644	Female / 14(Y)	Dr H
62546	TR20210507-0-001-0000		Sahiwal Lab_Main	2021-05-07	Iqra Shokat	6546564	Female / 30(Y)-9(M)-9(D)	Dr S
62522	TR20210506-0-004-0000		CC2	2021-05-06	Irfan Ahmad	0302562289	Male / 25(Y)-10(D)	Dr F
62521	TR20210506-0-001-0001		Sahiwal Lab_Main	2021-05-06	Ayza Khan	0302562289	Female / 22(Y)-3(M)-9(D)	Dr F
62520	TR20210506-0-001-0000		Sahiwal Lab_Main	2021-05-06	Khalid Ayub	654654	Male / 22(Y)-3(M)-24(D)	Dr F
62505	TR20210505-0-003-0001		CC1	2021-05-05	Hira Khan	654654654	Female / 22(Y)-3(M)-2(D)	
62502	TR20210505-0-004-0000		CC2	2021-05-05	Malik Shanzain	6565654	Male / 22(Y)-4(M)-7(D)	
62501	TR20210505-0-003-0000		CC1	2021-05-05	Hira Khan	654654654	Female / 22(Y)-3(M)-2(D)	Dr F
62500	TR20210505-0-001-0000		Sahiwal Lab_Main	2021-05-05	Khalid Ayub	654654	Male / 22(Y)-3(M)-24(D)	Dr F
62493	TR20210504-0-001-0002		Sahiwal Lab_Main	2021-05-04	Irfan Ahmad	0302562289	Male / 25(Y)-10(D)	Dr F
62492	TR20210504-0-003-0003		CC1	2021-05-04	Hira Khan	654654654	Female / 22(Y)-3(M)-2(D)	
62491	TR20210504-0-003-0002		CC1	2021-05-04	Ali Ahmad	554654	Male /	Dr S
62490	TR20210504-0-002-0003		Sahiwal_Collection_Ce...	2021-05-04	Khalid Ayub	654654	Male / 22(Y)-3(M)-24(D)	Dr L
62489	TR20210504-0-5586-0...		Testcc	2021-05-04	3rd Patient	65656	Male / 1(D)	
62488	TR20210504-0-5586-0...		Testcc	2021-05-04	Ayza Khan	0302562289	Female / 22(Y)-3(M)-9(D)	
62487	TR20210504-0-002-0002		Sahiwal_Collection_Ce...	2021-05-04	Ayza Khan	0302562289	Female / 22(Y)-3(M)-9(D)	
62484	TR20210504-0-003-0001		CC1	2021-05-04	Khalid Ayub	654654	Male / 22(Y)-3(M)-24(D)	Dr F
62483	TR20210504-0-003-0000		CC1	2021-05-04	Malik Shanzain	6565654	Male / 25(Y)-4(M)-7(D)	Dr S
62482	TR20210504-0-001-0001		Sahiwal Lab_Main	2021-05-04	Ali Ahmad	554654	Male / 22(Y)-3(M)-24(D)	Dr F
62481	TR20210504-0-004-0000		CC2	2021-05-04	Hira Khan	654654654	Female / 22(Y)-3(M)-2(D)	Dr F
62472	TR20210504-0-002-0001		Sahiwal_Collection_Ce...	2021-05-04	Patient 4	+5+65+65	Female / 30(Y)	

Sample settings

On samples listing a link displaying on the top left of the page on click a screen will display admin can set column which wants.

The screenshot shows the 'Sample listing setting' dialog box. It has a title bar with 'Sample listing setting' and a close button. The dialog contains a table with three columns: 'Column name', 'Display order', and 'Visible'. The 'Visible' column contains checkboxes. A 'Samples settings' link is highlighted in the background.

Column name	Display order	Visible
T/R Id	1	<input checked="" type="checkbox"/>
T/R No.	2	<input checked="" type="checkbox"/>
M/R No.	3	<input checked="" type="checkbox"/>
Collection centre	4	<input checked="" type="checkbox"/>
Collection date	5	<input checked="" type="checkbox"/>
Patient	6	<input checked="" type="checkbox"/>
Patient phone#	7	<input checked="" type="checkbox"/>
Patient type / Age	8	<input checked="" type="checkbox"/>
Doctor	9	<input checked="" type="checkbox"/>
Doctor share	10	<input checked="" type="checkbox"/>
Total amount	11	<input checked="" type="checkbox"/>
Discount(%)	12	<input checked="" type="checkbox"/>
Paid amount	13	<input checked="" type="checkbox"/>
Due amount	14	<input checked="" type="checkbox"/>
Created by	15	<input checked="" type="checkbox"/>
Modified date	16	<input checked="" type="checkbox"/>
Modified by	17	<input checked="" type="checkbox"/>
Lab	18	<input checked="" type="checkbox"/>

Admin should change the display order.

Test result

Add test result

Once tests have been conducted and results are ready, they can be added into the system using the 'Add Test Result' screen.

Use the 'Test reference number' field to search the test against which the results need to be added.

Open the result in edit mode and write the test result with appropriate details.

Use the 'Print result' button to print the results. Please note that results are only printed if a pathologist has already approved the results and the patient has cleared the dues.

Patient

Resend an email to the patient for a reminder if dues are not clear, or results are ready to collect from the Lab.

Doctor

Resend email to doctor with a PDF file of the patient's report.

Multiple prints click to print with header

Multiple prints on letterhead report will print without header

The screenshot displays a web interface for managing test results. At the top, there is a table with columns: Parameter, Result, Result in words, Abnormal, Male, Female, Child, Child2, Infant, New born, Min max value in words, Unit, and Parameter type. Below the table, a red message states: "Test is approved by : Dr. Farah Naz Tahir Amin".

Below the message are three sections, each with a text input field and a clear button (X):

- Notes:** A text input field with a clear button.
- Comments:** A text input field with a clear button.
- Method:** A text input field with a clear button.

At the bottom of the interface, there are five blue buttons:

- Print result
- Resend email to patient
- Resend email to doctor
- Multiple prints
- Multiple prints on letter head

Test results reports:

Test result report now admin can remove or add header and footer using toggle button.

Print
 Toggle header
 Toggle footer




PAF_Child_Lab

Cannt
Phone: 65654 65465 , Whatsapp: 65654

PATIENT REPORT

Patient No:  40897 Patient Name: Hamza Ali Age/Gender: 21 (Y) - 8 (M) - 13 (D) / Male Phone: 923440041445 CNIC:	Order No:  63343 Location: PAF_Child_Lab_Main Collection Date/Time: 25-Jun-2021 12:28 PM Specimen: Specimen collection type
--	--

Bacteriology	Result			
<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%; border: 1px solid black;">Test</td> <td style="width: 40%; border: 1px solid black;">Normal range</td> <td style="width: 30%; border: 1px solid black;">Unit</td> </tr> </table>	Test	Normal range	Unit	 63343 25-Jun-2021 12:49 PM
Test	Normal range	Unit		

Electronically verified report.No signature(s) required (Not valid for court)Please contact the lab in case of emergency. You can view your reports online at <http://192.168.4.119>

LAB • HORMONES • PCR • DIGITAL X-RAY • COLOUR DOPPLER ULTRASOUND • ECHOCARDIOGRAPHY • ECG

Printing Date/Time: 25-Jun-2021 12:47 PM
Generated By: PAF_Child_Lab PAF_Child_Lab
Activate Windows
Go to Settings to activate Windows.

If you want to remove the header, check the toggle buttons.

Print
 Toggle header
 Toggle footer



PAF_Child_Lab

Cannt
Phone: 65654 65465 , Whatsapp: 65654

PATIENT REPORT

Patient No:  40897 Patient Name: Hamza Ali Age/Gender: 21 (Y) - 8 (M) - 13 (D) / Male Phone: 923440041445 CNIC:	Order No:  63343 Location: PAF_Child_Lab_Main Collection Date/Time: 25-Jun-2021 12:28 PM Specimen: Specimen collection type
--	--

Bacteriology	Result			
<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%; border: 1px solid black;">Test</td> <td style="width: 40%; border: 1px solid black;">Normal range</td> <td style="width: 30%; border: 1px solid black;">Unit</td> </tr> </table>	Test	Normal range	Unit	 63343 25-Jun-2021 12:49 PM
Test	Normal range	Unit		

Electronically verified report.No signature(s) required (Not valid for court)Please contact the lab in case of emergency. You can view your reports online at <http://192.168.4.119>

LAB • HORMONES • PCR • DIGITAL X-RAY • COLOUR DOPPLER ULTRASOUND • ECHOCARDIOGRAPHY • ECG

Printing Date/Time: 25-Jun-2021 12:47 PM
Generated By: PAF_Child_Lab PAF_Child_Lab
Activate Windows
Go to Settings to activate Windows.

Test results

All sample collections which are added in OneLABz can be viewed on this screen.

ONE LABZ Dashboard Sample collection Test result Patient Accounting Stock Quality control Reports Settings Access control Admin Sahiwal Lab Sahiwal Lab									
Test result Test result settings									
Export									
Copy CSV Excel PDF Print									
Search by <--Please select a value--> Search value Search value Search Reset									
T/R Id.	Test result id	T/R No.	M/R No.	Test name	Collection centre	Patient	Patient phone#	Patient t	
62552	73940	TR20210517-0-001-0000		Bacter	Sahiwal Lab_Main	Muneeba Kafayat	9254644	Female / 1-	
62552	73941	TR20210517-0-001-0000		allergy	Sahiwal Lab_Main	Muneeba Kafayat	9254644	Female / 1-	
62546	73931	TR20210507-0-001-0000		allergy	Sahiwal Lab_Main	Iqra Shokat	6546564	Female / 3I	
62522	73887	TR20210506-0-004-0000		CBC	CC2	Irfan Ahmad	0302562289	Male / 25Y	
62521	73886	TR20210506-0-001-0001		allergy	Sahiwal Lab_Main	Ayza Khan	0302562289	Female / 2I	
62520	73884	TR20210506-0-001-0000		allergy	Sahiwal Lab_Main	Khalid Ayub	654654	Male / 22Y	
62520	73885	TR20210506-0-001-0000		Bacter	Sahiwal Lab_Main	Khalid Ayub	654654	Male / 22Y	
62505	73862	TR20210505-0-003-0001		CBC	CC1	Hira Khan	654654654	Female / 2I	
62502	73859	TR20210505-0-004-0000		Bacter	CC2	Malik Shahzain	6565654	Male / 25Y	
62501	73858	TR20210505-0-003-0000		allergy	CC1	Hira Khan	654654654	Female / 2I	
62500	73857	TR20210505-0-001-0000		Bacter	Sahiwal Lab_Main	Khalid Ayub	654654	Male / 22Y	
62493	73848	TR20210504-0-001-0002		allergy	Sahiwal Lab_Main	Irfan Ahmad	0302562289	Male / 25Y	
62493	73849	TR20210504-0-001-0002		Bacter	Sahiwal Lab_Main	Irfan Ahmad	0302562289	Male / 25Y	
62493	73850	TR20210504-0-001-0002		Genetics	Sahiwal Lab_Main	Irfan Ahmad	0302562289	Male / 25Y	
62492	73847	TR20210504-0-003-0003		allergy	CC1	Hira Khan	654654654	Female / 2I	
62491	73846	TR20210504-0-003-0002		allergy	CC1	Ali Ahmad	554654	Male /	
62490	73845	TR20210504-0-002-0003		CBC	Sahiwal_Collection_Ce...	Khalid Ayub	654654	Male / 22Y	
62489	73844	TR20210504-0-5586-0...		Genetics	Testcc	3rd Patient	69656	Male / 1(D)	
62488	73843	TR20210504-0-5586-0...		Genetics	Testcc	Ayza Khan	0302562289	Female / 2I	
62487	73842	TR20210504-0-002-0002		Genetics	Sahiwal_Collection_Ce...	Ayza Khan	0302562289	Female / 2I	
62484	73839	TR20210504-0-003-0001		allergy	CC1	Khalid Ayub	654654	Male / 22Y	

Pending perform

List of tests which are not performed yet display on this screen. Admin can copy, download CSV, PDF file can print on clicking buttons that are on the left side of the page. Search filters can apply that are on the right side of the page.

ONE LABZ Dashboard Sample collection Test result Patient Accounting Stock Quality control Reports Settings Access control Admin PAF_Child_Lab PAF_Chil...									
Test result Test result settings									
All test results Pending perform Pending approval									
Export									
Copy CSV Excel PDF Print									
Search by <--Please select a value--> Search value Search value Search Reset									
T/R Id.	Test result id	T/R No.	M/R No.	Test name	Collection centre	Patient	Patient		
63343	74924	TR20210625-0-001-0004		Histo	PAF_Child_Lab_Main	Hamza Ali	9234400-		
63343	74923	TR20210625-0-001-0004		Hematology	PAF_Child_Lab_Main	Hamza Ali	9234400-		
63343	74922	TR20210625-0-001-0004		Endocrinology	PAF_Child_Lab_Main	Hamza Ali	9234400-		
63343	74921	TR20210625-0-001-0004		DNA	PAF_Child_Lab_Main	Hamza Ali	9234400-		
63343	74920	TR20210625-0-001-0004		CH patho	PAF_Child_Lab_Main	Hamza Ali	9234400-		
63343	74919	TR20210625-0-001-0004		CBC	PAF_Child_Lab_Main	Hamza Ali	9234400-		
63343	74918	TR20210625-0-001-0004		blood test	PAF_Child_Lab_Main	Hamza Ali	9234400-		
63338	74912	TR20210625-0-001-0000		Bacteriology	PAF_Child_Lab_Main	Hamza Ali	9234400-		

Showing 25 records per page. 1-8 Of 8

<<Start 1 End>>

Pending approval

List of tests which are not performed yet display on this screen.

Admin can copy, download CSV, PDF file can print on clicking buttons that are on the left side of the page. Search filters can apply that are on the right side of the page.

Send WhatsApp:

On click this button can send WhatsApp to inform you about the result report.

Send sms:

If you want to send sms on your mobile phone can send on click this button.

The screenshot shows the CNELABZ test result form. At the top, there is a navigation bar with various menu items. Below the navigation bar, the form contains several input fields for parameters, results, and notes. At the bottom of the form, there are several buttons: 'Save', 'Save and print result', 'Print single result', 'Multiple prints', 'Send whatsapp', and 'Send sms'. The 'Send sms' button is highlighted with a red rectangle. There are also checkboxes for 'Show test history on report', 'Show header on report', and 'Show footer on report'. At the bottom right, there are buttons for 'Resend email to doctor' and 'Resend email to patient'. The footer of the page includes 'Copyright ©2021 Powered by TechXect' and a Windows activation watermark.

The screenshot shows the CNELABZ test result listing page. At the top, there is a navigation bar with various menu items. Below the navigation bar, the page title is 'Test result'. There are tabs for 'All test results', 'Pending perform', and 'Pending approval'. Below the tabs, there are export options (Copy, CSV, Excel, PDF, Print) and a search bar. The main content is a table of test results with columns: T/R Id., Test result id, T/R No., M/R No., Test name, Collection centre, Patient, and Patient phone#. The table contains 10 rows of data. At the bottom of the table, there is a 'Showing 25 records per page. 1-10 Of 10' indicator. At the bottom right, there is a 'Test result settings' link highlighted with a red rectangle. The footer of the page includes 'Activate Windows' and a Windows activation watermark.

Test result settings


On test result listing a link is displayed on the top of the page when click on it a page displaying admin wants to change column orders or no of grid use this.

The screenshot shows the CNELABZ test result listing settings page. The page title is 'Test result listing setting'. It contains a table with columns: Column name, Display order, and Visible. The table lists 13 columns and their corresponding display orders. The 'Test result settings' link from the previous page is highlighted with a red rectangle. The footer of the page includes 'Sahiwal Lab Sahiwal Lab' and a Windows activation watermark.

Column name	Display order	Visible
T/R Id.	1	<input checked="" type="checkbox"/>
Test result id	2	<input checked="" type="checkbox"/>
T/R No.	3	<input checked="" type="checkbox"/>
M/R No.	4	<input checked="" type="checkbox"/>
Test name	5	<input checked="" type="checkbox"/>
Collection centre	6	<input checked="" type="checkbox"/>
Patient	7	<input checked="" type="checkbox"/>
Patient phone#	8	<input checked="" type="checkbox"/>
Patient type / Age	9	<input checked="" type="checkbox"/>
Doctor	10	<input checked="" type="checkbox"/>
Collection date	11	<input checked="" type="checkbox"/>
Performed on	12	<input checked="" type="checkbox"/>
Approval status	13	<input checked="" type="checkbox"/>

Expense

Expense is below the Accounting menu.




Add Expense

There are three types of expenses managed in the system.

- Petty cash
- Purchases
- Salary


To add expenses, use the 'Expense' menu and select the 'Add expenses' option. Fill the form with relevant details and save the form to record an expense.



Expenses can be added for any collection center.

Expenses

All expenses which are added in OneLABz can be viewed on this screen.




The screenshot shows the 'Expenses' page in the OneLABz system. At the top, there is a navigation bar with the OneLABz logo and various menu items: Dashboard, Sample collection, Test result, Patient, Accounting, Stock, Quality control, Reports, Settings, Access control, and Admin. The user's name 'Inova lab Inova lab' is visible in the top right corner. Below the navigation bar, the page title 'Expenses' is displayed. There are several action buttons: 'Export' (with sub-options: Copy, CSV, Excel, PDF, Print), 'Search by' (a dropdown menu), 'Search value' (a text input field), and 'Search' and 'Reset' buttons. The main content is a table with the following columns: Expense id, Expense type, Amount (Rs.), Description, Expense date, Created date, Created by, Modified date, Modified by, Collection centre, and Lab. The table contains 11 rows of data. At the bottom of the table, there is a 'Showing 25 records per page, 1-10 Of 10' indicator and a pagination control with '«Start', '<', '1', '>', and 'End»' buttons.

Expense id	Expense type	Amount (Rs.)	Description	Expense date	Created date	Created by	Modified date	Modified by	Collection centre	Lab
30148	Petty Cash	1500.000	Fifteen hundred r...	19-04-2021	19-04-2021	Inova lab Inova lab	19-04-2021	Inova lab Inova lab	Inova lab_Main	Inova lab
20134	Salary	10000.000	Ten thousand only	24-03-2021	24-03-2021	Inova lab Inova lab	24-03-2021	Inova lab Inova lab	Chughtai	Inova lab
20133	Petty Cash	500.000	Five hundred only.	24-03-2021	24-03-2021	Inova lab Inova lab	24-03-2021	Inova lab Inova lab	Hum CC	Inova lab
20127	Petty Cash	5000.000	five thousand rup...	09-02-2021	18-03-2021	Inova lab Inova lab	18-03-2021	Inova lab Inova lab	Inova lab_Main	Inova lab
20118	Purchases	2000.000	two thousand only	15-03-2021	15-03-2021	Inova lab Inova lab	15-03-2021	Inova lab Inova lab	Jinnah road	Inova lab
20115	Petty Cash	1500.000	Fifteen hundred r...	15-03-2021	15-03-2021	Inova lab Inova lab	15-03-2021	Inova lab Inova lab	Inova lab_Main	Inova lab
20114	Purchases	300.000	Three hundred ru...	12-03-2021	12-03-2021	Inova Lab operator	12-03-2021	Inova Lab operator	Inova lab_Main	Inova lab
20113	Purchases	800.000	Eight hundred ru...	11-03-2021	12-03-2021	Inova lab Inova lab	12-03-2021	Inova lab Inova lab	Inova lab_Main	Inova lab
20112	Purchases	1000.000	One thousand ru...	12-03-2021	12-03-2021	Inova lab Inova lab	12-03-2021	Inova lab Inova lab	Inova lab_Main	Inova lab
20099	Petty Cash	10.000	ten rupees only	04-03-2021	04-03-2021	Inova lab Inova lab	04-03-2021	Inova lab Inova lab	Inova lab_Main	Inova lab

Expense reports:

There are four expense report




The screenshot shows the 'Expense report' form in the OneLABz system. The form has a navigation bar at the top with the OneLABz logo and menu items: Dashboard, Sample collection, Test result, Patient, Accounting, Stock, Quality control, Reports, Settings, Access control, and Admin. The user's name 'Inova lab Inova lab' is visible in the top right corner. The form contains several input fields and dropdown menus: 'Lab *' (dropdown menu with 'Inova lab' selected), 'Collection centre' (dropdown menu with 'Collection centre' selected), 'From' (text input field with 'DD-MM-YYYY' placeholder), 'To' (text input field with 'DD-MM-YYYY' placeholder), 'Order by field *' (dropdown menu with 'Date' selected), 'Order by' (dropdown menu with 'Ascending' selected), 'Report *' (dropdown menu with 'Daily expenses summary' selected), and 'Employee' (text input field with 'Employee' placeholder). There are three buttons at the bottom: 'View report', 'Print report', and 'Reset'.

You can view

- Daily expenses report
- Monthly summary report
- Profit and loss
- daily expenses details


Add expense type:

Add the expense type, fill all mandatory fields and click on the save button to save it.



Expense type:

All expense types which are added in OneLABz can be viewed on this screen.




Reports

There are many reports available in OneLABz for analysis and strategic decisions based on different data. To generate a report, use the 'Report' menu and select the particular report for which data is needed.

Patient dues amount report

This report is to figure out the due amount of patients.



There are certain filters available which can be used to filter report data

Patient dues report

Patient name: Aamir Sohail - 03076610559
Phone number: 03076610559
From: DD/MM/YYYY
To: DD/MM/YYYY

Show only unpaid tests:
Show due summary:

View report | Print report | Reset

Patient dues report

Total amount 0.000 | Total discount 0.000 | Total referral amount 0.000 | Total paid 0.000 | Total due 0.000

Doctor shares Report

Doctor shares report

Doctor name: Doctor name
Phone number: Doctor phone number
From: DD/MM/YYYY
To: DD/MM/YYYY

View report | Print report | Reset

Dropdown menu options:
Patient due amount report
Doctor shares report
Daily summary report
Doctor share summary report

A laboratory may agree referral amounts with doctors, for patients that the doctor refers to the lab.

There are certain filters available which can be used to filter report data.

ONELABZ Dashboard Sample collection Test result Patient Accounting Stock Quality control Reports Settings Access control Admin Sahiwal Lab Sahiwal Lab

Doctor shares report

Doctor name: From: To:

Order by field: Order by:

Doctor id

Doctor id	Doctor	Total share
20262	Dr Arshad	1092.000
20264	Dr Rahil	3611.440
20275	Dr Shahid	1164.864
20276	Dr Laila shfiq	97.300
20277	Dr Fg	730.000
20284	Dr Fazal	308.000
20302	Dr Shabir	120.000
20303	Dr Moeen	115.500

Copyright ©2021 Powered by TechXect

Daily summary report

The summary report provides a daily business summary. It includes Earnings, Expenses, Discounts, Paid amount, Balance due, Referral amount and Net amount for each date.

ONELABZ Dashboard Sample collection Test result Patient Accounting Stock Quality control Reports Settings Access control Admin Sahiwal Lab Sahiwal Lab

Daily summary report

From: To: Employee:

Collection centre:

Order by field: Order by:

View report Print report Reset

Date	Earnings	Expenses	Discounts	Paid amount	Previous amount	Balance due	Referral amount	Net amount
2021-04-17	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000
2021-04-18	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000
2021-04-19	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000
2021-04-20	3800.000	0.000	0.000	300.000	0.000	3800.000	280.000	-280.000
2021-04-21	36400.000	1500.000	2042.800	23457.200	0.000	5900.000	2796.440	24160.760
2021-04-22	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000
2021-04-23	0.000	0.000	0.000	5000.000	5000.000	0.000	0.000	0.000
2021-04-24	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000
2021-04-25	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000
2021-04-26	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000
2021-04-27	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.000
2021-04-28	2070.000	0.000	0.000	2070.000	0.000	0.000	0.000	0.000
2021-04-29	5000.000	0.000	0.000	5000.000	0.000	0.000	500.000	4500.000
2021-04-30	6250.000	0.000	0.000	6250.000	0.000	0.000	0.000	6250.000

Total amounts of all aspects are at the bottom of the page.

12/01/2020	20350.000	0.000	3660.055	16340.000	350.000	690.000	15649.945						
11/01/2020	23850.000	2810.000	2169.675	19730.000	1950.010	120.000	16800.315						
10/01/2020	6150.000	0.000	50.000	6100.000	0.000	0.000	6100.000						
09/01/2020	24750.000	0.000	4859.845	19140.000	749.970	180.000	18960.185						
08/01/2020	22250.000	0.000	889.990	21010.090	350.000	870.000	20140.010						
07/01/2020	14600.000	0.000	1590.075	9440.000	3570.000	540.000	8899.925						
06/01/2020	3050.000	0.000	50.000	2350.000	650.000	0.000	2350.000						
05/01/2020	15450.000	600.000	611.500	12298.500	2540.000	854.400	10844.100						
04/01/2020	20700.000	0.000	2390.000	12110.000	6200.000	600.000	11510.000						
03/01/2020	2950.000	0.000	0.000	2950.000	0.000	0.000	2950.000						
02/01/2020	9750.000	1000.000	106.500	9543.500	100.000	0.000	8543.500						
01/01/2020	6050.000	0.000	30.000	5970.000	50.000	0.000	5970.000						
31/12/2019	3550.000	0.000	600.000	2950.000	0.000	0.000	2950.000						
30/12/2019	650.000	0.000	0.000	650.000	0.000	0.000	650.000						
29/12/2019	6600.000	0.000	0.000	6600.000	0.000	0.000	6600.000						
26/12/2019	5300.000	0.000	0.000	5300.000	0.000	0.000	5300.000						
25/12/2019	700.000	0.000	0.000	0.000	700.000	0.000	0.000						
13/12/2019	150.000	0.000	0.000	150.000	0.000	0.000	150.000						
12/12/2019	750.000	0.000	100.000	650.000	0.000	195.000	455.000						
07/12/2019	650.000	0.000	0.000	650.000	0.000	0.000	650.000						
05/12/2019	3900.000	0.000	2250.000	0.000	1650.000	495.000	-495.000						
Earnings	332,500.000	Expenses	87,210.000	Discounts	31,857.611	Paid amount	260,332.100	Balance due	40,310.010	Referral amount	9,764.400	Net amount	163,357.979

Doctor Share Summary Report

This report shows a summary of a doctor's share for all tests performed on patients that were referred by doctors. Records where no doctor is specified are for patients that were not referred by a doctor with referral fee.

Doctor share summary report

Doctor name:

Patient name:

Doctor phone number:

Patient phone number:

Test name:

Lab id:

Date from:

Date to:

View report | Print report | Reset

Doctor share summary report

Lab id	Collection centre	Collection date	Doctor name	Patient name	Doctor share	Gross charges	Discount	Test Name
50935	Main lab	01/22/2020		Shagufta	0.0000	300.0000	0.0000	Semen Analysis Report
50934	Main lab	01/22/2020		Zafar Iqbal	0.0000	100.0000	0.0000	Blood Grouping

At the bottom of that page, total Gross charges, discounts, doctor's share are available.

30003	Main lab	12/26/2019		Wajid Ali	0.0000	300.0000	0.0000	Semen Analysis Report
20004	Main lab	12/25/2019	Dr.Farah Javed sb.	Maimoona Bibi	0.0000	400.0000	0.0000	CBC (Complete Blood Count)
20003	Main lab	12/25/2019		Mohammad Akram	0.0000	300.0000	0.0000	RFT (Renal Function Test)
10007	Main lab	12/13/2019		Shabbir Hussain	0.0000	50.0000	0.0000	Blood Glucose Report Fasting
10006	Main lab	12/13/2019		Sarwar Bibi	0.0000	100.0000	0.0000	Blood Glucose Report Fasting,Haemoglobin Level
10005	Main lab	12/12/2019		Mrs. Hafiz Abdul Hafeez	0.0000	50.0000	0.0000	Blood Glucose Report Fasting
10004	Main lab	12/12/2019		Mrs. Hafiz Abdul Hafeez	0.0000	50.0000	0.0000	
10003	Main lab	12/12/2019	Dr. Muhammad Rauf ch.	Ayesha Bibi	195.0000	650.0000	0.0000	CBC (Complete Blood Count),ESR (Erythrocyte Sedimentation Rate),Urine Complete Examination
10002	Main lab	12/07/2019	Dr. Muhammad Rauf ch.	Ali Ahmad	0.0000	650.0000	0.0000	CBC (Complete Blood Count),ESR (Erythrocyte Sedimentation Rate),Urine Complete Examination
3	Main lab	12/05/2019	Dr. Muhammad Rauf ch.	Shaisat Bibi	300.0000	1000.0000	0.0000	CBC (Complete Blood Count),H.Pylori Antibodies(Screening)
2	Main lab	12/05/2019	Dr. Muhammad Rauf ch.	Huma Bibi	195.0000	650.0000	0.0000	CBC (Complete Blood Count),ESR (Erythrocyte Sedimentation Rate),Urine Complete Examination
1	Main lab	12/05/2019		Mrs. Hafiz Abdul Hafeez	0.0000	2250.0000	0.0000	Lipid Profile,Blood Glucose Report Fasting,Uric Acid,LFT (Liver Function Test),RFT (Renal Function Test),Urine Complete Examination
Gross charges		330,000.000	Total discount	0.000	Doctor share	9,764,400		

Test Summary Report:

In this report detail of tests against lab and patient recorded.Admin can view patient and test detail using this report.

ONE LABS Dashboard Sample collection Test result Patient Accounting Stock Quality control Reports Settings Access control Admin Sahiwal Lab Sahiwal Lab

Test summary report

From: YYYY-MM-DD To: YYYY-MM-DD Employee: Employee

Test category: Test category Test: Test

Order by field *
 T/R No. (selected)
 <- Please select a value ->
 T/R No.
 Date
 Employee
 Patient name
 Tests category
 Test
 Earnings
 Discounts
 Balance due
 Paid amount
 Net amount

Order by: Ascending

Test summary report

T/R No.	From	To	Employee	Tests category	Test	Earnings	Discounts	Balance due	Paid amount	Net amount
6227				Bacteriology Blood Test	Bacter CBC	1000.0000	0.0000	1000.0000	0.0000	0.0000
62279	2021-04-20		Sahiwal Lab	Bacteriology	Bacter	500.0000	0.0000	500.0000	0.0000	0.0000
62283	2021-04-20		Sahiwal Lab	Bacteriology	Bacter	500.0000	0.0000	500.0000	0.0000	0.0000
62285	2021-04-20		Sahiwal Lab	Bacteriology Blood Test	Bacter CBC blood test	1300.0000	0.0000	1300.0000	0.0000	0.0000
62290	2021-04-20		Sahiwal Lab	Bacteriology	Bacter	500.0000	0.0000	500.0000	0.0000	0.0000

Payment Details Report:

In this report against TR/No Earnings, Balance due, Paid amount, Remaining amount, Doc share, payment date, Created date, Patient name, Test category, Test detail is present.

It can search TR/No.

Arrange in different order by field.

ONE LABS Dashboard Sample collection Test result Patient Accounting Stock Quality control Reports Settings Access control Admin Sahiwal Lab Sahiwal Lab

Payment detail report

From: YYYY-MM-DD To: YYYY-MM-DD Sample received by: Sample received by

Test category: Test category Test: Test Payment received by: Payment received by

T/R No. (input field)
 Order by field *
 T/R No. (selected)
 <- Please select a value ->
 T/R No.
 Earnings
 Balance due
 Paid amount
 Remaining amount
 Discounts
 Doc share
 Payment date
 Sample create date
 Patient name
 Payment received by
 Sample received by
 Test category
 Test

Order by: Ascending

View report Print report Reset

T/R No.	Earnings	Balance due	Paid amount	Remaining amount	Discounts	Doc share	Payment date	Sample create date	Patient name	Payment received by	Test category	Test
62278	1000.0000	1000.0000	0.0000	1000.0000	0.0000	0.0000			Patient_1_sahiwal Lab	Sahiwal Lab	Bacteriology Blood Test	Bacter CBC
62279	500.0000	500.0000	0.0000	500.0000	0.0000	0.0000	2021-04-20	2021-04-20	Patinet_3_sahiwal_Lab	Sahiwal Lab	Bacteriology	Bacter
62279	500.0000	500.0000	0.0000	500.0000	0.0000	0.0000	2021-04-20	2021-04-20	Patinet_3_sahiwal_Lab	Sahiwal Lab	Bacteriology	Bacter
62283	500.0000	500.0000	0.0000	500.0000	0.0000	0.0000	2021-04-20	2021-04-20	Patinet_3_sahiwal_Lab	Sahiwal Lab	Bacteriology	Bacter
62285	1300.0000	1300.0000	0.0000	1300.0000	0.0000	0.0000	2021-04-20	2021-04-20	3rd Patient	Sahiwal Lab	Bacteriology Blood Test	Bacter CBC blood test

Company rates report:

In this report Rates are defined according to the collection center. Admin can view the Commission, Test charges, Beared discount amount, payable.

ONE LABS Dashboard Sample collection Test result Patient Accounting Stock Quality control Reports Settings Access control Admin Sahiwal Lab Sahiwal Lab

Company rates report

From: YYYY-MM-DD To: YYYY-MM-DD Collection centre: Collection centre

Order by field: Test name Order by: Ascending

View report Print report Reset

Company rates report

Test name	Date	Commission	Test charges	Beared discount amount	Payable
allergy	2021-05-04	0.000	3600.000	0.000	3600.000
allergy	2021-05-05	0.000	600.000	0.000	600.000
allergy	2021-05-06	0.000	1200.000	0.000	1200.000
allergy	2021-05-07	0.000	600.000	0.000	600.000
allergy	2021-05-17	0.000	600.000	0.000	600.000
Bacter	2021-04-20	0.000	3500.000	0.000	3500.000
Bacter	2021-04-21	0.000	5500.000	0.000	5500.000
Bacter	2021-04-28	0.000	500.000	0.000	500.000
Bacter	2021-04-30	0.000	500.000	0.000	500.000
Bacter	2021-05-04	0.000	500.000	0.000	500.000
Bacter	2021-05-05	0.000	1000.000	0.000	1000.000
Bacter	2021-05-06	0.000	500.000	0.000	500.000
Bacter	2021-05-17	0.000	500.000	0.000	500.000
blood test	2021-04-20	0.000	810.000	0.000	510.000

Payment received amount:
Detail of payment against date and employee.

ONE LABS Dashboard Sample collection Test result Patient Accounting Stock Quality control Reports Settings Access control Admin Sahiwal Lab Sahiwal Lab

Payment received report

From: YYYY-MM-DD To: YYYY-MM-DD Employee: Employee

Order by field: Date Order by: Ascending

Payment received details report

Date	Received amount	Employee name	Expenses
2021-04-20	300.0000	Sahiwal Lab Sahiwal Lab	0.000
2021-04-21	23457.2000	Sahiwal Lab Sahiwal Lab	1500.000
2021-04-23	5000.0000	Sahiwal Lab Sahiwal Lab	0.000
2021-04-28	2070.0000	Sahiwal Lab Sahiwal Lab	0.000
2021-04-29	5000.0000	Sahiwal Lab Sahiwal Lab	0.000
2021-04-30	6250.0000	Sahiwal Lab Sahiwal Lab	0.000
2021-05-03	250.0000	Sahiwal Lab Sahiwal Lab	0.000
2021-05-04	17700.0000	Sahiwal Lab Sahiwal Lab	0.000
2021-05-05	1050.0000	Sahiwal Lab Sahiwal Lab	0.000
2021-05-06	2200.0000	Sahiwal Lab Sahiwal Lab	0.000
2021-05-07	600.0000	Sahiwal Lab Sahiwal Lab	0.000
2021-05-17	770.0000	Sahiwal Lab Sahiwal Lab	0.000

Stock consumption report:
By clicking the View button can view the Detail of stock consumption with chemical name, Quantity and Date.

ONE LABS Dashboard Sample collection Test result Patient Accounting Stock Quality control Reports Settings Access control Admin Sahiwal Lab Sahiwal Lab

Stock consumption report

From: YYYY-MM-DD To: YYYY-MM-DD Chemical / Item: Chemical / Item

Order by field *
 Chemical/Item
 <- Please select a value ->
 Chemical/Item
 Stock quantity consumed
 Date

Order by: Ascending

Stock consumption report

Chemical	Stock quantity consumed	Date
Blood Testing kit	36.0000	2021-04-20
Blood Testing kit	105.0000	2021-04-21
Blood Testing kit	6.0000	2021-04-28
Blood Testing kit	10.0000	2021-04-29
Blood Testing kit	16.0000	2021-04-30
Blood Testing kit	76.0000	2021-05-04
Blood Testing kit	10.0000	2021-05-05
Blood Testing kit	5.0000	2021-05-06
Blood Testing kit	5.0000	2021-05-17
Chlorine	60.0000	2021-05-04
Chlorine	10.0000	2021-05-05
Chlorine	20.0000	2021-05-06
Chlorine	10.0000	2021-05-07
Chlorine	10.0000	2021-05-17

Doctor business report:

Can view the details of doctors' referral tests and the amount of every test its total.

Metro City Lab Dashboard Sample collection Test result Patient Accounting Stock Quality control Reports Settings Access control Admin Metro City Lab

Doctor business report

From: DD/MM/YYYY To: DD/MM/YYYY Employee * Metro City Lab - 03 400 18 5555

Test category: Test category Doctor: Doctor Collection centre: Collection centre

Order by field * Bacteriology Order by: Ascending

View report Print report Reset

Doctor business report

Employee: Metro City Lab - 03 400 18 5555

Collection centre:

DoctorName	Blood Test	DNA	Sugar	Microbiology	Clinical Chemistry	Heamatology	Molecular Biology	Histopathology	Immunology	Bacteriology	Endocrinology	Toxicology	Serology	Chemical Pathology	XRAY	MRI	None	Special Chemistry	SC-1	Total amount
10 March Patient	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000
11 March Doctor	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000
14-April-SMS	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000

Discount option report

In the discount option report detail of discount is present. The user can view the date, discount option, employee and collection center.

Discount option report

From: To: Discount option:

Employee: Collection centre:

Order by field *: Order by:

Discount option report

DiscountOptionName	Gross amount	Discount amount	Paid amount	Due amount
Discount 2	3450.0000	690.0000	2760.0000	0.0000
Discount 3	46000.0000	18400.0000	27600.0000	0.0000
Discount for Regular customer	10770.0000	3231.0000	7539.0000	0.0000